

# AMPLATZER™ VASCULAR PLUGS

Embolization Therapy

PRECISE PLACEMENT AND RAPID  
EMBOLIZATION WITH A SINGLE DEVICE<sup>1,2,3</sup>



Amplatzer™  
Vascular Plug (AVP)



Amplatzer™  
Vascular Plug II (AVP II)



Amplatzer™  
Vascular Plug 4 (AVP 4)

1. Pech M, Kraetsch A, Wieners G, et al. Embolization of the Gastroduodenal Artery Before Selective Internal Radiotherapy: A prospectively Randomized Trial Comparing Platinum-Fibered Microcoils with the AMPLATZER Vascular Plug II. CVIR 2009(32)3:455-61.

2. Kucukay F, Özdemir M, Şenol E, Okten S, Ereren M, Karan A. Large pulmonary arteriovenous malformations: long-term results of embolization with AMPLATZER vascular plugs. J Vasc Interv Radiol. 2014 Sep;25(9):1327-32. doi: 10.1016/j.jvir.2014.01.031. Epub 2014 Mar 18.

3. Tests performed by and data on file at Abbott.

**Information contained herein for DISTRIBUTION outside the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

©2020 Abbott. All rights reserved. MAT-2007616 v1.0



# AMPLATZER™ FAMILY OF VASCULAR PLUGS

## Peripheral Vascular Embolization

### ADVANCING THE STANDARD OF CARE IN PERIPHERAL EMBOLIZATION WITH ABBOTT

- A single device solution that allows for faster procedure times, less radiation exposure and lower procedural costs<sup>1,2</sup>
- Rapid embolization even in high-flow vessels<sup>1,3,4</sup>
- Designed for precise placement, controlled delivery and fully recapturable
- A family of plugs for different vessel conditions



1. Pech M, Kraetsch A, Wieners G, et al. Embolization of the Gastroduodenal Artery Before Selective Internal Radiotherapy: A prospectively Randomized Trial Comparing Platinum-Fibered Microcoils with the AMPLATZER Vascular Plug II. CVIR 2009(32)3:455-61.

2. Jackson J, Hart J, Aldin Z, et al. Embolization of pulmonary arteriovenous malformations using the AMPLATZER vascular plug: successful treatment of 69 consecutive patients. Eur Radiol 2010;20(11)2663-70.

3. Kucukay F, Özdemir M, Şenol E, Okten S, Ereren M, Karan A. Large pulmonary arteriovenous malformations: long-term results of embolization with AMPLATZER vascular plugs. J Vasc Interv Radiol. 2014 Sep;25(9):1327-32. doi: 10.1016/j.jvir.2014.01.031. Epub 2014 Mar 18.

4. Tests performed by and data on file at Abbott.

**Information contained herein for DISTRIBUTION outside the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.**




# AMPLATZER™ VASCULAR PLUG

## Short Landing Zone Embolization

**Compact design:** Single-lobe nitinol mesh design ideal for short landing zones<sup>1</sup>


**Guide catheter or sheath deliverable:** Compatible with 4-6F sheaths or 5-8F guide catheters depending on device size



Treatable vessel size range:  
2.7 - 12.3 mm<sup>a</sup>

## SIZING AND DEVICE SELECTION

Vessel size	Device dimensions			Delivery catheter requirements <sup>2</sup>				
Treatable vessel diameter range <sup>†</sup>	Model / reorder number	Vascular plug diameter [A]	Unconstrained length [B]	Minimum internal diameter	Minimum sheath Size	or	Minimum guide catheter size	Maximum delivery catheter length <sup>3</sup>
2.5 mm - 3.0 mm	9-PLUG-004	4 mm	7 mm	≥ 1.42 mm / ≥ 0.056 in	≥ 4 F		≥ 5 F	≤ 100 cm
4.0 mm - 4.5 mm	9-PLUG-006	6 mm	7 mm	≥ 1.42 mm / ≥ 0.056 in	≥ 4 F		≥ 5 F	≤ 100 cm
5.5 mm - 6.0 mm	9-PLUG-008	8 mm	7 mm	≥ 1.42 mm / ≥ 0.056 in	≥ 4 F		≥ 5 F	≤ 100 cm
6.5 mm - 7.5 mm	9-PLUG-010	10 mm	7 mm	≥ 1.68 mm / ≥ 0.066 in	≥ 5 F		≥ 6 F	≤ 100 cm
8.0 mm - 9.0 mm	9-PLUG-012	12 mm	8 mm	≥ 1.68 mm / ≥ 0.066 in	≥ 5 F		≥ 6 F	≤ 100 cm
9.5 mm - 11.0 mm	9-PLUG-014	14 mm	8 mm	≥ 2.21 mm / ≥ 0.087 in	≥ 6 F		≥ 8 F	≤ 100 cm
10.5 mm - 12.5 mm	9-PLUG-016	16 mm	8 mm	≥ 2.21 mm / ≥ 0.087 in	≥ 6 F		≥ 8 F	≤ 100 cm



<sup>†</sup> Treatable vessel diameter range based on the devices Instructions for Use to select a plug that is oversized by approximately 30-50% at the occlusion site.

<sup>a</sup> Treatable vessel size range based on the device's Instructions for Use to select a device with a diameter approximately 30-50% larger than the vessel diameter at the occlusion site.

1. Lopera, Jorge E. "The Amplatzer vascular plug: review of evolution and current applications." Seminars in interventional radiology. Vol. 32. No. 04. Thieme Medical Publishers, 2015.

2. The AMPLATZER Vascular Plug is delivered utilizing either a sheath or guide catheter meeting the minimum internal diameter requirements above.

3. Each AMPLATZER Vascular Plug comes pre-loaded on a 135cm nitinol delivery wire.

**Information contained herein for DISTRIBUTION outside the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.**




# AMPLATZER™ VASCULAR PLUG II

## Variable Landing Zone Embolization<sup>1</sup>

**Rapid embolization<sup>2</sup>:** Multi-layered, multiple-lobed nitinol mesh design provides for rapid embolization within the vessel

**Secure positioning:** Multiple points of contact with the vessel wall for secure positioning in medium- and high-flow vessels


**Guide catheter or sheath deliverable<sup>3</sup>:** Compatible with 4-7F sheaths or 5-9F guide catheters depending on device size



Treatable vessel size range:  
2.0 - 16.9 mm<sup>a</sup>

## SIZING AND DEVICE SELECTION

Vessel size	Device dimensions			Delivery catheter requirements			
Treatable vessel diameter range†	Model / reorder number	Vascular plug diameter [A]	Unconstrained length [B]	Minimum and maximum internal diameter	Minimum sheath size	or <sup>3</sup> Minimum guide catheter Size	Maximum delivery catheter length <sup>4</sup>
2.0 mm - 2.5 mm	9-AVP2-003	3 mm	6 mm	1.42 - 1.70 mm / 0.056 - 0.067 in	≥ 4 F	≥ 5 F	≤ 100 cm
2.5 mm - 3.0 mm	9-AVP2-004	4 mm	6 mm	1.42 - 2.50 mm / 0.056 - 0.098 in	≥ 4 F	≥ 5 F	≤ 100 cm
4.0 mm - 4.5 mm	9-AVP2-006	6 mm	6 mm	1.42 - 2.50 mm / 0.056 - 0.098 in	≥ 4 F	≥ 5 F	≤ 100 cm
5.5 mm - 6.0 mm	9-AVP2-008	8 mm	7 mm	1.42 - 2.69 mm / 0.056 - 0.106 in	≥ 4 F	≥ 5 F	≤ 100 cm
6.5 mm - 7.5 mm	9-AVP2-010	10 mm	7 mm	1.78 - 2.69 mm / 0.070 - 0.106 in	≥ 5 F	≥ 6 F	≤ 100 cm
8.0 mm - 9.0 mm	9-AVP2-012	12 mm	9 mm	1.78 - 2.69 mm / 0.070 - 0.106 in	≥ 5 F	≥ 6 F	≤ 100 cm
9.5 mm - 11.0 mm	9-AVP2-014	14 mm	10 mm	2.18 - 2.69 mm / 0.086 - 0.106 in	≥ 6 F	≥ 8 F	≤ 100 cm
10.5 mm - 12.5 mm	9-AVP2-016	16 mm	12 mm	2.18 - 2.69 mm / 0.086 - 0.106 in	≥ 6 F	≥ 8 F	≤ 100 cm
12.0 mm - 14.0 mm	9-AVP2-018	18 mm	14 mm	2.49 - 2.69 mm / 0.098 - 0.106 in	≥ 7 F	≥ 9 F	≤ 100 cm
13.5 mm - 15.5 mm	9-AVP2-020	20 mm	16 mm	2.49 - 2.69 mm / 0.098 - 0.106 in	≥ 7 F	≥ 9 F	≤ 100 cm
14.5 mm - 17.0 mm	9-AVP2-022	22 mm	18 mm	2.49 - 2.69 mm / 0.098 - 0.106 in	≥ 7 F	≥ 9 F	≤ 100 cm



a. Treatable vessel size range based on the device's Instructions for Use to select a device with a diameter approximately 30-50% larger than the vessel diameter at the occlusion site. † Treatable vessel diameter range based on the devices Instructions for Use to select a plug that is oversized by approximately 30-50% at the occlusion site. 1. AVP II is available in 8 unconstrained lengths ranging from 6 -18 mm. 2. Pech M., Kraetsch A., Wieners G., et al. Embolization of the Gastroduodenal Artery Before Selective Internal Radiotherapy: A Prospectively Randomized Trial Comparing Platinum-Fibred Microcoils with the AMPLATZER Vascular Plug II. Cardiovasc Intervent Radiol. 2009;32(3)455-61. 3. The AMPLATZER Vascular Plug II is delivered utilizing either a sheath or guide catheter meeting the minimum internal diameter requirements above. 4. Each AMPLATZER Vascular Plug II comes pre-loaded on a 135cm nitinol delivery wire.

**Information contained herein for DISTRIBUTION outside the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.**



# AMPLATZER™ VASCULAR PLUG 4


## Low-Profile Embolization

**Diagnostic catheter deliverable:** Simple delivery through a 0.038” guide wire-compatible diagnostic catheter

**Improved navigation:** Low-profile design and more flexible delivery wire allow the device to navigate through tortuous anatomies with ease<sup>2</sup>


**Rapid embolization<sup>1</sup>:** Multi-layered, double-lobed nitinol mesh design provides for rapid embolization within the vessel

Treatable vessel size range: 2.6 - 6.2 mm<sup>a</sup>



## SIZING AND DEVICE SELECTION

Vessel size	Device dimensions			Diagnostic catheter requirements	
Treatable vessel diameter range <sup>†</sup>	Model / reorder number	Vascular plug diameter [A]	Unconstrained length [B]	Requirements for diagnostic catheters <sup>3,4</sup>	Maximum delivery catheter length <sup>5</sup>
2.5 - 3.0 mm	9-AVP038-004	4 mm	10.0 mm	0.038” Guide Wire-Compatible Diagnostic Catheter	≤ 125 cm
3.5 - 4.0 mm	9-AVP038-005	5 mm	10.5 mm	0.038” Guide Wire-Compatible Diagnostic Catheter	≤ 125 cm
4.0 - 4.5 mm	9-AVP038-006	6 mm	11.0 mm	0.038” Guide Wire-Compatible Diagnostic Catheter	≤ 125 cm
4.5 - 5.5 mm	9-AVP038-007	7 mm	12.5 mm	0.038” Guide Wire-Compatible Diagnostic Catheter	≤ 125 cm
5.5 - 6.0 mm	9-AVP038-008	8 mm	13.5 mm	0.038” Guide Wire-Compatible Diagnostic Catheter	≤ 125 cm



<sup>†</sup> Treatable vessel diameter range based on the devices Instructions for Use to select a plug that is oversized by approximately 30-50% at the occlusion site.

<sup>a</sup> Treatable vessel size range based on the device's Instructions for Use to select a device with a diameter approximately 30-50% larger than the vessel diameter at the occlusion site.

<sup>1</sup> Tests performed by and data on file at Abbott.

<sup>2</sup> Lopera, Jorge E. “The Amplatzer vascular plug: review of evolution and current applications.” Seminars in interventional radiology. Vol. 32. No. 04. Thieme Medical Publishers, 2015.

<sup>3</sup> The AMPLATZER Vascular Plug 4 is delivered utilizing an 0.038” Guidewire-Compatible Diagnostic Catheter with adequate wall strength.

<sup>4</sup> The AMPLATZER Vascular Plug 4 has been tested for compatibility with the following diagnostic catheters and corresponding lengths:

5F Dx Catheters: Boston Scientific IMAGER II (≤ 100 cm in length) and Merit Medical Impress (≤ 125 cm in length)

4F Dx Catheters: Cordis TEMPO or TEMPO AQUA (≤ 100 cm in length)

<sup>5</sup> Each AMPLATZER Vascular Plug 4 comes pre-loaded and packaged on an 0.038” x 155cm-long PTFE-coated delivery wire.

**Information contained herein for DISTRIBUTION outside the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.**



The AMPLATZER™ Vascular Plug, AMPLATZER™ Vascular Plug II, and AMPLATZER™ Vascular Plug 4 are each indicated for arterial and venous embolizations in the peripheral vasculature.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at [www.vascular.eifu.abbott](http://www.vascular.eifu.abbott) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

**Information contained herein for DISTRIBUTION outside the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

**Abbott International BVBA**  
Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

[www.cardiovascular.abbott](http://www.cardiovascular.abbott)  
©2020 Abbott. All rights reserved. MAT-2007616 v1.0

