Read Highlighted Changes: Revised November 2015.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

ARCHITECT Anti-CCP

INTENDED USE

The ARCHITECT Anti-CCP assay is a chemiluminescent microparticle immunoassay (CMIA) for the semi-quantitative determination of the IgG class of autoantibodies specific to cyclic citrullinated peptide (CCP) in human serum or plasma on the ARCHITECT iSystem. Detection of anti-CCP antibodies is used as an aid in the diagnosis of Rheumatoid Arthritis (RA) and should be used in conjunction with other clinical information. Autoantibody levels represent one parameter in a multicriterion diagnostic process, encompassing both clinical and laboratory-based assessments.

SUMMARY AND EXPLANATION OF THE TEST

Rheumatoid Arthritis (RA) is a common, systemic autoimmune disease affecting 0.5-1% of the population. It is characterized by chronic inflammation of the synovium, which commonly leads to progressive joint destruction and in most cases, to disability and reduction of quality of life.¹ Evidence gained over the last few years suggests that aggressive therapy given early in the disease has the greatest therapeutic potential.^{2, 3}

The serum of RA patients contains a variety of antibodies directed against self-antigens. The most widely known of these autoantibodies is the rheumatoid factor (RF) antibody directed against the constant domain of IgG molecules. The presence of RF is one of the American College of Rheumatology's (ACR) criteria for the classification of RA.⁴ Although the RF test has good sensitivity for RA, it is not very specific for the disease as it can also be detected in the serum of patients with other rheumatic or inflammatory diseases and even in a substantial percentage of the healthy (elderly) population.⁵ For several years it has been recognized that antibodies to anti-perinuclear factor (APF) and anti-keratin (AKA) are highly specific for RA. It was subsequently reported that both of these antibodies reacted with native filaggrin and are now referred to as anti-filaggrin antibodies (AFA).⁶⁻⁸ More recently it has been shown that all of these antibodies are directed to citrullinecontaining epitopes.⁹ Citrulline is a non-standard amino acid, as it is not incorporated into proteins during protein synthesis. It can, however, be generated via post-translational modification of arginine residues by the enzyme peptidyl arginine deiminase (PAD).¹⁰ In 1998, Schellekens and colleagues reported that linear peptides containing citrulline (CP) were very specific for RA antibodies (96%) in an ELISA based assay.¹¹ Subsequent work demonstrated that cyclic variants of these peptides, termed cyclic citrullinated peptides (CCP), were equally specific for RA, but with a higher sensitivity than linear peptides.¹² To improve the sensitivity of the CCP test further, several dedicated libraries of citrulline-containing peptides were screened with RA sera and a new set of peptides (CCP2) were discovered which gave superior performance compared to the CCP1 test.¹³ Over the last few years, many independent studies have confirmed the diagnostic performance of the CCP2 test.^{14, 15} In 2007, the European League against Rheumatism (EULAR) published guidelines for the diagnosis of early RA, and the measurement of antibodies to anti-CCP was included as a serology marker.¹⁶

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT Anti-CCP assay is a two-step immunoassay with an automated sample pretreatment for the semi-quantitative determination of the IgG class of autoantibodies specific to cyclic citrullinated peptide (CCP) in human serum or plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

- Sample is prediluted with wash buffer. The prediluted sample, CCP coated paramagnetic microparticles and sample diluent are combined. The Anti-CCP antibodies present in the sample binds to the CCP coated microparticles.
- 2. After washing, anti-human IgG acridinium-labeled conjugate is added to create a reaction mixture.
- 3. Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
- The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of anti-CCP antibody in the sample and the RLUs detected by the ARCHITECT iSystem optics.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Kit Contents

ARCHITECT Anti-CCP 1P65

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

REF	1P65-25	1P65-35
Σ	100	500
MICROPARTICLES	1 x 6.5 mL	1 x 26.5 mL
CONJUGATE	1 x 5.8 mL	1 x 25.8 mL
SAMPLE DILUENT	1 x 9.8 mL	1 x 50.0 mL

MICROPARTICLES CCP coated microparticles in phosphate buffer with surfactant and protein (bovine) stabilizer. Minimum concentration: 0.05% solids. Preservative: sodium azide.

CONJUGATE Mouse anti-human IgG: acridinium-labeled conjugate in MES buffer with surfactant and protein (bovine) stabilizer. Minimum concentration: 10 ng/mL. Preservatives: Nipasept and Sarafloxacin.

SAMPLE DILUENT Phosphate buffer with surfactant and protein (bovine) stabilizer. Preservative: sodium azide.

Other Reagents

PRE-TRIGGER SOLUTION ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

TRIGGER SOLUTION ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

WASH BUFFER ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

Warnings and Precautions

- IVD
- For In Vitro Diagnostic Use

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁷⁻²⁰

The following warnings an	nd precautions apply to: MICROPARTICLES
and SAMPLE DILUENT	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in
	accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagent kits beyond the expiration date.
- Do not pool reagents within a kit or between kits.
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE**, Assay Procedure section of this package insert.
- Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.
 - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
 - Once a septum has been placed on an open reagent bottle, do not invert the bottle as this will result in reagent leakage and may compromise assay results.
 - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

When stored and handled as directed, reagents are stable until the expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened/ Opened*	2-8°C	Until expiration date	May be used immediately after removal from 2-8°C storage. Store in upright position.
On board	System temperature	30 days	Discard after 30 days. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the** microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded. For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The ARCHITECT Anti-CCP assay file must be installed on the ARCHITECT iSystem prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

Alternate Result Units

The default result unit for the ARCHITECT Anti-CCP assay is U/mL.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

Verified specimen types to be used with this assay:

Specimen Types	Collection Tubes	
Serum	Serum	
	Serum separator tubes	
Plasma	Lithium heparin plasma separator	
	tubes	
	Potassium EDTA	

- Other specimen collection tube types have not been tested with this assay.
- Plasma specimens from different anticoagulant tube types should not be used interchangeably for monitoring anti-CCP.
- Performance has not been established for the use of cadaveric specimens or the use of body fluids other than human serum or plasma.
- Liquid anticoagulants may have a dilution effect resulting in lower concentrations for individual patient specimens.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
 - pooled
 - grossly hemolyzed
 - obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

• Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.

- Mix thawed specimens thoroughly by low speed vortexing or by inverting 10 times. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous.
- To ensure consistency in results, specimens must be transferred to a centrifuge tube and centrifuged at ≥ 10,000 RCF (Relative Centrifugal Force) for 10 minutes before testing if
 - they contain fibrin, red blood cells, or other particulate matter,
 - they require repeat testing, or
 - they were frozen and thawed.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

•	-	
Specimen Type	Storage Temperature	Maximum Storage Time
Serum/Plasma	Room temperature (study performed at 30°C)	≤ 22 hours
	2-8°C	≤ 7 days
	-20°C or colder	

Specimens may be stored on or off the clot, red blood cells, or separator gel.

If testing will be delayed more than 22 hours for specimens stored at room temperature or more than 7 days for specimens stored at 2-8°C, remove serum or plasma from the clot, red blood cells, or separator gel and store at -20°C or colder.

Avoid more than three freeze/thaw cycles.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

PROCEDURE

Materials Provided

1P65 ARCHITECT Anti-CCP Reagent Kit

Materials Required but not Provided

- ARCHITECT Anti-CCP Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbottdiagnostics.com.
- 1P65-01 ARCHITECT Anti-CCP Calibrators
- 1P65-10 ARCHITECT Anti-CCP Controls
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
 - Invert the microparticle bottle 30 times.

- Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
- If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.
- Once the microparticles have been resuspended, place a septum on the bottle. For instructions about placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles.
- Order calibration, if necessary.
 - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
 - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

Maximum number of replicates sampled from the same sample cup: 10

• Priority:

Sample volume for first test: 60 µL

- Sample volume for each additional test from same sample cup: 10 μL
- \leq 3 hours on board:

Sample volume for first test: 150 µL

Sample volume for each additional test from same sample cup: 10 μL

- If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT Anti-CCP Calibrators and Controls.
 - ARCHITECT Anti-CCP Calibrators and Controls should be prepared according to their respective package inserts.
 - Recommended volumes:
 - for each calibrator: 4 drops
 - for each control: 4 drops
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Specimen Dilution Procedures

Specimens with anti-CCP value exceeding 200.0 U/mL are flagged with the code "> 200.0 U/mL" and may be diluted using either the Automated Dilution Protocol or the Manual Dilution Procedure.

Automated Dilution Protocol

The system performs a 1:6 dilution of the specimen and automatically calculates the concentration of the specimen before dilution and reports the result.

Manual Dilution Procedure

Suggested dilution: 1:10.

1. Add 50 μL of the patient specimen to 450 μL of ARCHITECT Anti-CCP Negative Control.

 The operator must enter the dilution factor in the Patient or Control order screen. The system will use this dilution factor to automatically calculate the concentration of the sample before dilution and report the result.

For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

Calibration

 Test Calibrators A-F in replicates of two. The calibrators should be priority loaded.

A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.

- Calibration Range: 0.0 200.0 U/mL.
- Once an ARCHITECT Anti-CCP calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

The recommended control requirement for the ARCHITECT Anti-CCP assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.

Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

Control values must be within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated test results are invalid and must be retested. Recalibration may be indicated.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

The ARCHITECT Anti-CCP assay belongs to method group 1.

RESULTS

Calculation

The ARCHITECT Anti-CCP assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration curve.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

Measurement Range (Reportable Range)

The measurement range of the ARCHITECT Anti-CCP assay is 0.5 U/mL to 200.0 U/mL.

LIMITATIONS OF THE PROCEDURE

- For diagnostic purposes, the ARCHITECT Anti-CCP results should be used in conjunction with other clinical data; e.g., symptoms, medical history, etc.
- If the ARCHITECT Anti-CCP results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- The value of anti-CCP in juvenile arthritis has not been determined.
- Some specimens may not dilute linearly because of the heterogeneity of the autoantibodies with respect to physiochemical properties.
- ARCHITECT Anti-CCP results should not be used interchangeably with other manufacturers' methods for anti-CCP determinations.

- Plasma specimens from different anticoagulant tube types should not be used interchangeably for monitoring anti-CCP.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA).^{21, 22} Specimens containing HAMA may produce anomalous values when tested with assay kits such as ARCHITECT Anti-CCP that employ mouse monoclonal antibodies.²¹
- Heterophilic antibodies in human specimens can react with reagent immunoglobulins, interfering with *in vitro* immunoassays.²³ Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed. Additional information may be required for diagnosis.
- Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS section in this package insert for specimen limitations.

EXPECTED VALUES

In a representative study, serum specimens from 199 asymptomatic, apparently healthy males (n=126) and females (n=73), with an age range of 19 to 67 years, were tested with the ARCHITECT Anti-CCP assay. No differences attributable to gender or age were observed. Specimen values ranged from < 0.5 U/mL to 2.5 U/mL. A cut-off of 5.0 U/mL was chosen, whereby a result of \geq 5.0 U/mL is considered positive and a result of <5.0 U/mL is considered negative.*

* Representative data; results in individual laboratories may vary from these data. It is recommended that each laboratory establish its own expected range, which may be unique to the population it serves depending upon geographical, patient, dietary, or environmental factors.

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The ARCHITECT Anti-CCP assay is designed to have an imprecision of < 10% total CV.

A study was performed based on guidance from the National Committee for Clinical Laboratory Standards (NCCLS) document EP5-A2.²⁴ Seven samples consisting of the ARCHITECT Anti-CCP Positive Control, four human plasma panels, and two human plasma samples were assayed on two instruments, in replicates of two at two separate times per day for 20 days (n = 80 for each sample), using two lots of reagents and a single calibration for each instrument/ reagent lot combination. Data from this study are summarized in the following table.*

		Reagent		Mean	Withi	n Run	То	tal
Sample	Instrument	Lot	n	(U/mL)	SD	%CV	SD	%CV
Positive	1	1	80	24.5	0.73	3.0	0.81	3.3
Control	2	2	80	26.7	0.68	2.6	0.74	2.8
Panel 1	1	1	80	10.9	0.30	2.7	0.61	5.6
	2	2	80	11.3	0.26	2.3	0.58	5.2
Panel 2	1	1	80	28.6	0.63	2.2	1.95	6.8
	2	2	80	30.3	0.60	2.0	1.68	5.5
Panel 3	1	1	80	66.7	1.40	2.1	4.03	6.0
	2	2	80	72.7	1.92	2.6	4.85	6.7
Panel 4	1	1	80	135.3	6.36	4.7	8.11	6.0
	2	2	80	154.1	5.02	3.3	11.82	7.7
Sample 1	1	1	80	2.8	0.07	2.6	0.11	4.0
	2	2	80	2.7	0.07	2.4	0.11	4.0
Sample 2	1	1	80	181.4	6.67	3.7	9.69	5.3
	2	2	80	195.3	7.20	3.7	12.14	6.2

* Representative data; results in individual laboratories may vary from these data.

Sensitivity

Sensitivity is defined as the limit of detection (LoD). The ARCHITECT Anti-CCP assay is designed to have a LoD of \leq 0.5 U/mL. The LoD and the limit of blank (LoB) of the ARCHITECT Anti-CCP assay were determined based on guidance from the NCCLS document EP17-A²⁵ using proportions of false positives (a) less than 5% and false negatives (β) less than 5%. These determinations were performed using one blank (60 replicates) and five low level anti-CCP samples (20 replicates each); LoB = 0.02 U/mL and LoD = 0.11 U/mL.*

* Representative data; results in individual laboratories may vary from these data.

Linearity

The ARCHITECT Anti-CCP assay is designed to be linear across the measurement range of 0.5 to 200.0 U/mL.

Based on a study performed by guidance from the NCCLS document EP6-A,^{26} the ARCHITECT Anti-CCP assay demonstrated linearity from 0.5 to 200.0 U/mL.*

* Representative data; results in individual laboratories may vary from these data.

Concentration Range (U/mL)	Slope (95% CI)	Intercept (95% CI)	r ²
0.1 - 257.4	0.98	-1.85	0.9985
	(0.95 to 1.01)	(-6.19 to 2.48)	

Autodilution Verification

The ARCHITECT Anti-CCP automated dilution method is designed to have a mean difference of $\pm 10\%$ versus the manual dilution method when performed on samples with values > 50.0 U/mL.

The ARCHITECT Anti-CCP assay was evaluated with the 1:6 autodilution method versus the 1:10 manual dilution method using 12 human serum samples with anti-CCP levels ranging from 58.7 to 785.0 U/mL. Five replicates each of the autodiluted and manually diluted samples were assayed on one instrument using the ARCHITECT Anti-CCP assay. The mean percent difference across all samples was 2.6%. The percent difference results are summarized in the following table.*

	Mean Automated Diluted Value x Dilution Factor of 6	Mean Manually Diluted Value x Dilution Factor of 10	
Sample	(U/mL)	(U/mL)	% Difference ^a
1	456.4	453.0	0.7
2	504.1	482.7	4.4
3	796.8	743.8	7.1
4	734.6	785.0	-6.4
5	220.2	209.9	4.9
6	192.0	187.9	2.2
7	213.9	207.0	3.3
8	196.8	194.6	1.1
9	65.3	61.4	6.3
10	70.1	69.2	1.3
11	72.0	69.0	4.4
12	167.2	165.2	1.2
	[Mean Auton	nated Diluted Value x 6 (U/mL) -	
a 0/ Difforon	Mean Manu	ally Diluted Value x 10 (U/mL)]	
/ DILIFIELD			X 100

Mean Manually Diluted Value x 10 (U/mL)

* Representative data; results in individual laboratories may vary from these data.

Interference

The ARCHITECT Anti-CCP assay is designed to have a maximum deviation in anti-CCP concentration from the following potentially interfering compounds within:

- ±15% for anti-CCP concentrations > 10.0 U/mL
- $\pm 10\%$ for anti-CCP concentrations ≥ 5.0 U/mL to ≤ 10.0 U/mL
- ±0.5 U/mL for anti-CCP concentrations < 5.0 U/mL

A study was performed based on guidance from the Clinical and Laboratory Standards Institute (CLSI) document EP7-A2²⁷ for the ARCHITECT Anti-CCP assay. Serum samples with anti-CCP levels across the assay range of 0.5 U/mL to 200.0 U/mL were supplemented with the potentially interfering compounds listed in the table below. The maximum deviation of anti-CCP concentration observed in serum samples during these studies ranged from:

- -7.6% to 0.8% for anti-CCP concentrations > 10.0 U/mL
- -1.0% to 7.5% for anti-CCP concentrations \geq 5.0 U/mL to \leq 10.0 U/mL
- -0.3 U/mL to 0.2 U/mL for anti-CCP concentrations < 5.0 U/mL*

	No Interference Found up to the Following
Potentially Interfering Substance	Concentration
Bilirubin	20 mg/dL
Hemoglobin	800 mg/dL
Total Protein	12 g/dL
Triglycerides	3000 mg/dL
Rheumatoid Factor	200 IU/mL
Red Blood Cells	0.4%

* Representative data; results in individual laboratories may vary from these data.

Cross-Reactivity

To assess the potential cross-reactivity of the CCP antigen used in the ARCHITECT Anti-CCP assay with other autoantibodies, the assay was evaluated with 20 samples positive for various other autoantibodies and negative for CCP antibodies. The following autoantibodies (1-4 samples of each) were tested in the assay: SSA, SSB, Sm, RNP, ds-DNA, Jo-1, ScI-70, Ribo-P, TPO, ANA, and AMA. The study showed no significant cross-reactivity of the CCP antigen with any of these other autoantibodies.

Tube Type Matrix Comparison

The specimen collection tubes listed below were verified for use with the ARCHITECT Anti-CCP assay.

serum, serum separator, lithium heparin plasma separator, and potassium EDTA.

When compared to the control tube type (serum), the tube types evaluated for samples with anti-CCP values < 5.0 U/mL showed less than a 0.5 U/mL difference on average, and the tube types evaluated for samples with anti-CCP values ranging from 5.3 to 178.8 U/mL showed less than a 10% difference on average. The distribution of the differences or percent differences per tube type is listed in the following table.*

	Distribution of Absolute Differences < 0.5 U/mL for Samples with Anti-	Distribution of Absolute Percent Differences for Samples with Anti-CO Values 5.3 to 178.8 U/mL		rcent Anti-CCP 1L	
Tube Type	CCP Values < 5.0 U/mL	$< 10\%$ $\geq 10\%$ to $\leq 20\%$ > 20\%			
Serum Separator	100%	76%	16%	8%	
	(19/19)	(19/25)	(4/25)	(2/25)	
Potassium EDTA	100%	72%	24%	4%	
	(19/19)	(18/25)	(6/25)	(1/25)	
Lithium Heparin	100%	80%	16%	4%	
Plasma Separator	(19/19)	(20/25)	(4/25)	(1/25)	

* Representative data; results in individual laboratories may vary from these data.

Clinical Sensitivity and Specificity

The clinical sensitivity was determined for 496 confirmed RA individuals, and clinical specificity was determined for 499 non-RA specimens (299 from patients with other rheumatic and non-rheumatic disorders and 200 from asymptomatic apparently healthy individuals). Using a cut-off of 5.0 U/mL, the sensitivity was calculated to be 70.6% with a specificity of 98.2%. The results are summarized in the following tables.*

	ARCHITECT Anti-CCP		
Specimen Category	Total n	Positive n	% Sensitivity
Confirmed RA ^a	496	350	70.6

^a RA patients were classified according to the ACR Criteria.⁴

	ARCHITECT Anti-CCP		
Specimen Category	Total n	Positive n	% Specificity
Non-RA Specimens in Total	499	9	98.2
Non-RA Healthy Asymptomatic	200	1	99.5
Non-RA Disease Specimens ^a	299	8	97.3

^a The non-RA diseases were Ankylosing Spondylitis, Autoimmune Thyroiditis/Hashimoto's Disease, Crohn's Disease, Dermatomyositis, Epstein-Barr Virus, Lyme Disease, Osteoarthritis, Polymyalgia Rheumatica, Polymyositis, Psoriatic Arthritis, Reactive Arthritis/ Reiter's Syndrome, Scleroderma, Sjögren's Syndrome, Systemic Lupus Erythematosus, and Ulcerative Colitis.

* Representative data; results in individual laboratories may vary from these data.

Method Comparison

The ARCHITECT Anti-CCP assay is designed to have a concordance of \geq 95% for RA and non-RA specimens when compared to the AxSYM Anti-CCP assay. The RA and non-RA specimens described in the Clinical Sensitivity and Specificity section were used to compare the ARCHITECT Anti-CCP assay to the AxSYM Anti-CCP assay. The cut-off employed for the AxSYM Anti-CCP assay was 5.0 U/mL, as stated in the manufacturer's package insert. Using a cut-off of 5.0 U/mL for the ARCHITECT Anti-CCP assay, the concordance was calculated to be 99.3%. The results are summarized in the following tables.*

			AxSYM Anti-CCP	
n = 995			Positive	Negative
ARCHITECT Anti-CCP		Positive	356	3
		Negative	4	632
	Total n	% Positive Agreement (95% Cl ^a)	% Negative Agreement (95% Cl ^a)	% Total Agreement (95% Cl ^a)
Non-RA	499	100 (66.4-100)	100 (99.2-100)	100 (99.3-100)
RA	496	98.9 (97.1-99.7)	97.9 (94.1-99.6)	98.6 (97.1-99.4)
All Samples	995	98.9 (97.2-99.7)	99.5 (98.6-99.9)	99.3 (98.6-99.7)

^a CI = Confidence Interval.

* Representative data; results in individual laboratories may vary from these data.

A Receiver Operator Characteristic (ROC) analysis was carried out using the above data obtained for the two assays. The area under the curve (AUC) for the ARCHITECT Anti-CCP assay was 0.873 (95% confidence interval: 0.849-0.897) and 0.872 (95% confidence interval: 0.848-0.896) for the AxSYM Anti-CCP assay, thus indicating that both assays are comparable with respect to their clinical differentiation. The ROC analysis curve is shown below.*



* Representative data; results in individual laboratories may vary from these data.

High Dose Hook

High dose hook is a phenomenon whereby very high level specimens may read within the dynamic range of the assay. For the ARCHITECT Anti-CCP assay, no high dose hook effect was observed when a sample containing approximately 2000 U/mL of anti-CCP antibody was assayed.*

* Representative data; results in individual laboratories may vary from these data.

BIBLIOGRAPHY

- Feldmann M, Brennan FM, Maini RN. Rheumatoid arthritis. Cell 1996;85:307-310
- Landewé RB. The benefits of early treatment in rheumatoid arthritis: confounding by indication, and the issue of timing. *Arthritis Rheum* 2003;48(1):1-5.
- Lard LR, Visser H, Speyer I, et al. Early versus delayed treatment in patients with recent-onset rheumatoid arthritis: comparison of two cohorts who received different treatment strategies. *Am J Med* 2001;111:446-451.
- Arnett FC, Edworthy SM, Bloch DA, et al. The American Rheumatism Association 1987 revised criteria for the classification of rheumatoid arthritis. *Arthritis Rheum* 1988;31(3):315-324.
- van Venrooij WJ, Hazes JM, Visser H. Anticitrullinated protein/ peptide antibody and its role in the diagnosis and prognosis of early rheumatoid arthritis. *Neth J Med* 2002;60(10):383-388.
- Nienhuis RL, Mandema E, Smids C. A new serum factor in patients with rheumatoid arthritis. The antiperinuclear factor. *Ann Rheum Dis* 1964;23:302-305.
- Young BJ, Mallya RK, Leslie RD, et al. Anti-keratin antibodies in rheumatoid arthritis. Br Med J 1979;2:97-99.
- Hoet RM, Boerbooms AM, Arends M, et al. Antiperinuclear factor, a marker autoantibody for rheumatoid arthritis: colocalisation of the perinuclear factor and profilaggrin. *Ann Rheum Dis* 1991;50:611-618.
- Sebbag M, Simon M, Vincent C, et al. The antiperinuclear factor and the so-called antikeratin antibodies are the same rheumatoid arthritis-specific autoantibodies. J Clin Invest 1995;95:2672-2679.
- Vossenaar ER, Zendman AJ, van Venrooij WJ, et al. PAD, a growing family of citrullinating enzymes: genes, features and involvement in disease. *BioEssays* 2003;25:1106-1118.
- Schellekens GA, de Jong BA, van den Hoogen FH, et al. Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies. *J Clin Invest* 1998;101(1):273-281.
- Schellekens GA, Visser H, de Jong BA, et al. The diagnostic properties of rheumatoid arthritis antibodies recognizing a cyclic citrullinated peptide. *Arthritis Rheum* 2000;43(1):155-163.
- Vossenaar ER, van Venrooij WJ. Anti-CCP antibodies, a highly specific marker for (early) rheumatoid arthritis. *Clin Applied Immunol Rev* 2004;4:239-262.
- Pruijn GJ, Vossenaar ER, Drijfhout JW, et al. Anti-CCP antibody detection facilitates early diagnosis and prognosis of rheumatoid arthritis. *Current Rheumatology Reviews* 2005;1(1):1-7.
- Nishimura K, Sugiyama D, Kogata Y, et al. Meta-analysis: diagnostic accuracy of anti-cyclic citrullinated peptide antibody and rheumatoid factor for rheumatoid arthritis. *Ann Intern Med* 2007;146(11):797-808.
- Combe B, Landewe R, Lukas C, et al. EULAR recommendations for the management of early arthritis: report of a task force of the European Standing Committee for International Clinical Studies Including Therapeutics (ESCISIT). Ann Rheum Dis 2007;66:34-45.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in* Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. Laboratory Biosafety Manual. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- 21. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.

- 22. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. Cancer Res 1985;45(2):879-885.
- 23. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34(1):27-33.
- 24. National Committee for Clinical Laboratory Standards (NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
- 25. National Committee for Clinical Laboratory Standards (NCCLS). Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline. NCCLS Document EP17-A. Wayne, PA: NCCLS; 2004.
- 26. National Committee for Clinical Laboratory Standards (NCCLS). Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. NCCLS Document EP6-A. Wayne, PA: NCCLS; 2003.
- 27. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI Document EP7-A2. Wayne, PA: CLSI; 2005.

Key to Symbols

i	Consult instructions for use	
	Manufacturer	
Σ	Sufficient for	
X	Temperature limitation	
	Use by/Expiration date	
CONJUGATE	Conjugate	
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.	
	Distributed in the USA by	
	Information needed for United	
INFORMATION FOR USA UNLI	States of America only	
IVD	<i>In Vitro</i> Diagnostic Medical Device	
LOT	Lot Number	
MICROPARTICLES	Microparticles	
PRE-TRIGGER SOLUTION	Pre-Trigger Solution	
PRODUCED FOR ABBOTT BY	Produced for Abbott by	
PRODUCT OF UK	Product of United Kingdom	
REACTION VESSELS	Reaction Vessels	
REAGENT LOT	Reagent Lot	
REF	List Number	
REPLACEMENT CAPS	Replacement Caps	
SAMPLE CUPS	Sample Cups	
SAMPLE DILUENT	Sample Diluent	
SEPTUM	Septum	
SN	Serial number	
TRIGGER SOLUTION	Trigger Solution	
WASH BUFFER	Wash Buffer	

ARCHITECT and AxSYM are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.



Abbott GmbH & Co. KG Max-Planck-Ring 2 65205 Wiesbaden Germany +49-6122-580

PRODUCED FOR ABBOTT BY

Abbott Laboratories

Axis-Shield Diagnostics Ltd.,

Dundee, UK

Abbott Park, IL 60064 USA

DISTRIBUTED IN THE USA BY

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised November 2015. ©2008, 2015 Abbott Laboratories



CE