



Read Highlighted Changes: Revised November 2016.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

**WARNING:** The Abbott ARCHITECT CA 19-9XR CMIA assay utilizes an antibody/antigen system based on the 1116-NS-19-9 antibody. The unique reagent formulation employed in the ARCHITECT CA 19-9XR assay may return elevated concentrations when compared to other methods for samples expressing high levels of 1116-NS-19-9 reactive determinants.<sup>1, 2</sup> Additionally, there is no internationally recognized standard for CA 19-9, which can contribute to differences between assay methods. The ARCHITECT CA 19-9XR assay is standardized to a reference standard prepared by Fujirebio Diagnostics, Inc. Performance characteristics of the Abbott ARCHITECT CA 19-9XR assay are NOT transferable to other diagnostic kits.

The concentration of 1116-NS-19-9 reactive determinants obtained with different assay methods cannot be used interchangeably due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the CA 19-9 assay used. If, in the course of monitoring a patient, the assay method used for determining serial 1116-NS-19-9 reactive determinant levels is changed, additional sequential testing should be carried out. Prior to changing assays, the laboratory MUST confirm baseline values for patients being serially monitored.

**WARNING:** 1116-NS-19-9 reactive determinants are shed naturally in saliva and other body fluids.<sup>3</sup> Contamination of the samples or the ARCHITECT iSystem disposables with saliva or aerosols (e.g., as a result of sneezing) may cause falsely elevated CA 19-9 assay values. It is recommended that all elevated values be reviewed and testing repeated as appropriate. Gloves should always be worn when handling samples, sample cups, reaction vessels, and septums. Face masks are also recommended.

## ■ NAME

ARCHITECT CA 19-9XR

## ■ INTENDED USE

The ARCHITECT CA 19-9XR assay is a chemiluminescent microparticle immunoassay (CMIA) for the quantitative determination of 1116-NS-19-9 reactive determinants in human serum or plasma on the ARCHITECT iSystem. The ARCHITECT CA 19-9XR assay is to be used as an aid in the management of pancreatic cancer patients in conjunction with other clinical methods.

## ■ SUMMARY AND EXPLANATION OF THE TEST

The ARCHITECT CA 19-9XR assay detects a tumor-associated antigen, which occurs in tissue as a monosialoganglioside and in serum as a high molecular weight, carbohydrate-rich glycoprotein known as a mucin.<sup>4-7</sup>

The ARCHITECT CA 19-9XR assay is based upon a monoclonal antibody, 1116-NS-19-9, which reacts with a carbohydrate antigenic determinant expressed on the circulating antigen.<sup>4-6</sup>

The results of published research studies<sup>8-14</sup> indicate that the CA 19-9 assay value is frequently elevated in the serum of subjects with various gastrointestinal conditions, such as pancreatic, colorectal, gastric, and hepatic carcinomas. No data exist to support the use of CA 19-9 in screening for malignancies.<sup>15, 16</sup> The role of CA 19-9 is to be used as an adjunct with other diagnostic information in the management of patients with pancreatic cancer.<sup>15</sup> Increased serum CA 19-9 assay values have also been observed in patients with metastases and in nonmalignant conditions such as hepatitis, cirrhosis, pancreatitis, and other gastrointestinal disease.<sup>8-11, 17-20</sup> Elevated levels have also been seen in cystic fibrosis.<sup>21-24</sup> Research studies demonstrate that CA 19-9 assay values may have utility in monitoring subjects with the above-mentioned diagnosed gastrointestinal malignancies.<sup>25-28</sup> It has been shown that a persistent elevation in CA 19-9 assay value following treatment may be indicative of occult metastatic and/or residual disease. A persistently rising CA 19-9 assay value may be associated with progressive malignant disease and poor therapeutic response. A declining CA 19-9 assay value may be indicative of a favorable prognosis and a good response to treatment.<sup>29-35</sup>

Testing for 1116-NS-19-9 reactive determinants must not be used as a screening procedure for malignancy. 1116-NS-19-9 reactive determinants are present as a normal constituent in serum and plasma of individuals without gastrointestinal carcinomas or having certain aforementioned non-cancer related conditions.

## ■ BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT CA 19-9XR assay is a two-step immunoassay for the quantitative determination of 1116-NS-19-9 reactive determinants in human serum or plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

1. Sample and 1116-NS-19-9 coated paramagnetic microparticles are combined. The 1116-NS-19-9 reactive determinants present in the sample bind to the 1116-NS-19-9 coated microparticles.
2. After washing, 1116-NS-19-9 acridinium-labeled conjugate is added to create a reaction mixture.
3. Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
4. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of 1116-NS-19-9 reactive determinants in the sample and the RLUs detected by the ARCHITECT iSystem optics.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

## REAGENTS

### Kit Contents

ARCHITECT CA 19-9XR 2K91

**NOTE:** Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

| REF                   | 2K91-32  | 2K91-24    | 2K91-39     |
|-----------------------|--|------------|-------------|
|                       | 100  | 400        | 500         |
| <b>MICROPARTICLES</b> | 1 x 6.6 mL   | 4 x 6.6 mL | 1 x 27.0 mL |
| <b>CONJUGATE</b>      | 1 x 5.9 mL   | 4 x 5.9 mL | 1 x 26.3 mL |
| <b>MICROPARTICLES</b> | 1116-NS-19-9 (mouse, monoclonal) coated microparticles in citrate buffer with protein (bovine) stabilizer. Minimum concentration: 0.09% solids. Preservatives: sodium azide and ProClin 300.       |            |             |
| <b>CONJUGATE</b>      | 1116-NS-19-9 (mouse, monoclonal) acridinium-labeled conjugate in phosphate buffer with protein (bovine) stabilizer. Minimum concentration: 0.5 µg/mL. Preservatives: sodium azide and ProClin 300. |            |             |

### Other Reagents

**MULTI-ASSAY MANUAL DILUENT** 1 x 100 mL ARCHITECT Multi-Assay Manual Diluent, **REF** 7D82-50, containing phosphate buffered saline solution. Preservative: antimicrobial agent.

**PRE-TRIGGER SOLUTION** ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

**TRIGGER SOLUTION** ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

**WASH BUFFER** ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

### Warnings and Precautions

- IVD**
- For *In Vitro* Diagnostic Use

### Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.<sup>36-39</sup>

| The following warnings and precautions apply to: <b>MICROPARTICLES</b> |  |
|--|--|
|  |  |
| <b>WARNING</b>   | Contains methylisothiazolones and sodium azide.                        |
| H317   | May cause an allergic skin reaction.                                   |
| EUH032   | Contact with acids liberates very toxic gas.                           |
| <b>Prevention</b>  |  |
| P261   | Avoid breathing mist / vapors / spray.                                 |
| P272   | Contaminated work clothing should not be allowed out of the workplace. |
| P280   | Wear protective gloves / protective clothing / eye protection.         |
| <b>Response</b>  |  |
| P302+P352  | IF ON SKIN: Wash with plenty of water.                                 |
| P333+P313  | If skin irritation or rash occurs: Get medical advice / attention.     |
| P362+P364  | Take off contaminated clothing and wash it before reuse.               |

| <b>Disposal</b>   |  |
|---|--|
| P501  | Dispose of contents / container in accordance with local regulations.  |
| The following warnings and precautions apply to: <b>CONJUGATE</b> |  |
|   |  |
| <b>DANGER</b>   | Contains polyethylene glycol octylphenyl ether, methylisothiazolones and sodium azide.   |
| H317  | May cause an allergic skin reaction.   |
| H318  | Causes serious eye damage.   |
| H412  | Harmful to aquatic life with long lasting effects.   |
| EUH032  | Contact with acids liberates very toxic gas.   |
| <b>Prevention</b>   |  |
| P261  | Avoid breathing mist / vapors / spray.   |
| P280  | Wear protective gloves / protective clothing / eye protection.   |
| P272  | Contaminated work clothing should not be allowed out of the workplace.   |
| P273  | Avoid release to the environment.  |
| <b>Response</b>   |  |
| P302+P352   | IF ON SKIN: Wash with plenty of water.   |
| P333+P313   | If skin irritation or rash occurs: Get medical advice / attention.   |
| P305+P351+P338  | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P310  | Immediately call a POISON CENTER or doctor / physician.  |
| P362+P364   | Take off contaminated clothing and wash it before reuse.   |
| <b>Disposal</b>   |  |
| P501  | Dispose of contents / container in accordance with local regulations.  |

Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

### Reagent Handling

- 1116-NS-19-9 reactive determinants are shed naturally in saliva and other body fluids.<sup>3</sup> Contamination of the samples or the ARCHITECT iSystem disposables with saliva or aerosols (e.g., as a result of sneezing) may cause falsely elevated CA 19-9 assay values. It is recommended that all elevated values be reviewed and testing repeated as appropriate. Gloves should always be worn when handling samples, sample cups, reaction vessels, and septums. Face masks are also recommended.**
- Do not use reagent kits beyond the expiration date.
- Do not pool reagents within a kit or between kits.**
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**

- To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
- Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
- Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

### Reagent Storage

When stored and handled as directed, reagents are stable until the expiration date.

|                         | Storage Temperature | Maximum Storage Time  | Additional Storage Instructions   |
|-------------------------|---------------------|-----------------------|---|
| <b>Unopened/Opened*</b> | 2-8°C               | Until expiration date | May be used immediately after removal from 2-8°C storage.<br>Store in upright position.   |
| <b>On board</b>         | System temperature  | 30 days               | Discard after 30 days.<br>For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5. |

\* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright.

**If any reagent bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

### Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

## INSTRUMENT PROCEDURE

The ARCHITECT CA 19-9XR assay file must be installed on the ARCHITECT iSystem prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

## SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

### Specimen Types

Verified specimen types to be used with this assay:

| Specimen Types | Collection Tubes      |
|----------------|-----------------------|
| Serum          | Serum                 |
|                | Serum separator tubes |
|                | Tripotassium EDTA     |
| Plasma         | Sodium Heparin        |
|                | Lithium Heparin       |

- Other specimen collection tube types have not been tested with this assay.

- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.
- Follow the tube manufacturer's processing instructions for serum and plasma collection tubes.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.

### Specimen Conditions

- Do not use specimens with the following conditions:
  - grossly hemolyzed
  - obvious microbial contamination
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.
- Performance has not been established using body fluids other than human serum or plasma.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

### Preparation for Analysis

- Follow the tube manufacturer's processing instructions for specimen collection tubes.
- Mix thawed specimens thoroughly by low speed vortexing or by inverting 10 times. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous.
- To ensure consistency in results, centrifuge specimens before testing if
  - they contain fibrin, red blood cells, or other particulate matter,
  - they require repeat testing, or
  - they were frozen and thawed.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

### Specimen Storage

| Specimen Type | Storage Temperature | Maximum Storage Time |
|---------------|---------------------|----------------------|
| Serum/Plasma  | 2-8°C               | ≤ 7 days             |
|               | -20°C or colder     | >7 days              |

- If testing will be delayed more than 24 hours, remove serum or plasma from the clot, serum separator or red blood cells.
- Specimens may be stored for up to 7 days at 2-8°C prior to being tested.
- If testing will be delayed more than 7 days, serum or plasma should be stored frozen at -20°C or colder.
- Avoid multiple freeze/thaw cycles.

### Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

## PROCEDURE

### Materials Provided

2K91 ARCHITECT CA 19-9XR Reagent Kit

### Materials Required but not Provided

- ARCHITECT CA 19-9XR Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com).
- 2K91 ARCHITECT CA 19-9XR Calibrators
- 2K91 ARCHITECT CA 19-9XR Controls
- 7D82-50 ARCHITECT Multi-Assay Manual Diluent
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

### Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
  - **Invert the microparticle bottle 30 times.**
  - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
  - **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
  - Once the microparticles have been resuspended, place a septum on the bottle. For instructions about placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
  - Verify that all necessary reagents are present.
  - Ensure that septums are present on all reagent bottles.
- Order calibration, if necessary.
  - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
  - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

Maximum number of replicates sampled from the same sample cup: 10

- Priority:
  - Sample volume for first test: 80  $\mu$ L
  - Sample volume for each additional test from same sample cup: 30  $\mu$ L
- $\leq$  3 hours on board:
  - Sample volume for first test: 150  $\mu$ L
  - Sample volume for each additional test from same sample cup: 30  $\mu$ L
- $>$  3 hours on board: Additional sample volume required

- If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT CA 19-9XR Calibrators and Controls.
  - Mix calibrator(s) and controls by gentle inversion before use.
  - Hold bottles **vertically** and dispense recommended volumes into each respective sample cup.
  - Recommended volumes:
    - for each calibrator: 4 drops
    - for each control: 4 drops
- Load samples.
  - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

### Specimen Dilution Procedures

Specimens with an ARCHITECT CA 19-9XR value exceeding 1200 U/mL are flagged with the code "> 1200.00" and may be diluted using either the Automated Dilution Protocol or the Manual Dilution Procedure.

#### Automated Dilution Protocol

The system performs a 1:10 dilution of the specimen and automatically calculates the concentration of the specimen before dilution and reports the result.

#### Manual Dilution Procedure

Suggested dilution: 1:10

An additional 1:10 dilution may be made if needed.

1. Add 50  $\mu$ L of the patient specimen to 450  $\mu$ L of ARCHITECT Multi-Assay Manual Diluent (7D82).
2. The operator must enter the dilution factor in the Patient or Control order screen. The system will use this dilution factor to automatically calculate the concentration of the sample before dilution and report the result. The dilution should be performed so that the diluted result reads greater than 30 U/mL.

For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

### Calibration

- Test Calibrators A-F in duplicate. The calibrators should be priority loaded.
  - A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Calibration Range: 0 - 1200 U/mL.
- Once an ARCHITECT CA 19-9XR calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
  - A reagent kit with a new lot number is used or
  - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

### Quality Control Procedures

The recommended control requirement for the ARCHITECT CA 19-9XR assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.

The ARCHITECT CA 19-9XR values must be within the acceptable ranges specified in the control package insert. If a control is out of the specified range, the associated test results are invalid and samples must be retested. Recalibration may be indicated.

## Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

The ARCHITECT CA 19-9XR assay belongs to method group 1.

## RESULTS

### Calculation

The ARCHITECT CA 19-9XR assay utilizes a Linear Regression data reduction method to generate a calibration curve.

### Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

## LIMITATIONS OF THE PROCEDURE

- The ARCHITECT CA 19-9XR assay value must be used in conjunction with information available from clinical evaluation and other diagnostic procedures.
- If the ARCHITECT CA 19-9XR results are inconsistent with clinical evidence, additional testing is recommended.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.<sup>40</sup>
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Specimens containing HAMA may produce anomalous values when tested with assay kits that employ mouse monoclonal antibodies. ARCHITECT CA 19-9XR reagents contain a component that reduces the effect of HAMA reactive specimens. Additional clinical or diagnostic information may be required to determine patient status.<sup>41, 42</sup>
- Patients with confirmed carcinoma may have pretreatment CA 19-9 assay values in the same range as healthy individuals. Elevations in circulating 1116-NS-19-9 reactive determinants may be observed in patients with metastases and in nonmalignant conditions such as hepatitis, cirrhosis, pancreatitis, and other gastrointestinal disease. Elevated levels have also been seen in cystic fibrosis.<sup>21</sup> For these reasons, a CA 19-9 assay value, regardless of level, should not be interpreted as absolute evidence for the presence or absence of malignant disease. **The ARCHITECT CA 19-9XR assay must not be used as a cancer screening test.**
- Patients with the Le<sup>a-b-</sup> phenotype may not express the 1116-NS-19-9 reactive determinant.<sup>43</sup>
- Representative performance data are given in the **EXPECTED VALUES and SPECIFIC PERFORMANCE CHARACTERISTICS** sections. Results obtained in individual laboratories may vary.

## EXPECTED VALUES

### APPARENTLY HEALTHY SUBJECTS

A study was performed with three hundred sixty (360) serum specimens from apparently healthy individuals. The distribution of ARCHITECT CA 19-9XR assay values from these specimens is shown in the table below.\*

|                             | Distribution of ARCHITECT CA 19-9XR Values |               |                 |                  |                   |              |
|-----------------------------|--|---------------|-----------------|------------------|-------------------|--------------|
|                             | Number of Subjects                         | Percent (%)   |                 |                  |                   |              |
|                             |  | 0-37.0 (U/mL) | 37.1-100 (U/mL) | 100.1-500 (U/mL) | 500.1-1200 (U/mL) | >1200 (U/mL) |
| Apparently Healthy Subjects | 360  | 94.4          | 5.6             | 0.0              | 0.0               | 0.0          |

In this study, 94.4% of the specimens from apparently healthy subjects (n=360) had values of 37 U/mL or less.

\* Representative data; results in individual laboratories may vary from these data.

## NONMALIGNANT DISEASE

A study was performed with four hundred forty one (441) samples from patients with nonmalignant disease to determine the distribution of serum ARCHITECT CA 19-9XR assay values. The distribution of values determined in this study is shown in the table below.\*

|                        | Number of Subjects | Distribution of ARCHITECT CA 19-9XR Values |                 |                  |                   |              |
|------------------------|--------------------|--|-----------------|------------------|-------------------|--------------|
|                        |                    | Percent (%)                                |                 |                  |                   |              |
| Nonmalignant Disease   |                    | 0-37.0 (U/mL)                              | 37.1-100 (U/mL) | 100.1-500 (U/mL) | 500.1-1200 (U/mL) | >1200 (U/mL) |
| Rectal Polyps          | 33                 | 97.0                                       | 3.0             | 0.0              | 0.0               | 0.0          |
| Pancreatitis           | 3                  | 100.0                                      | 0.0             | 0.0              | 0.0               | 0.0          |
| Gallbladder            | 21                 | 95.2                                       | 0.0             | 0.0              | 0.0               | 4.8          |
| Diabetes               | 38                 | 94.7                                       | 5.3             | 0.0              | 0.0               | 0.0          |
| Pulmonary              | 40                 | 100.0                                      | 0.0             | 0.0              | 0.0               | 0.0          |
| Cirrhosis              | 153                | 92.8                                       | 4.6             | 0.7              | 0.7               | 1.3          |
| Hepatitis              | 68                 | 92.6                                       | 7.4             | 0.0              | 0.0               | 0.0          |
| Renal                  | 34                 | 91.2                                       | 8.8             | 0.0              | 0.0               | 0.0          |
| Other Gastrointestinal | 51                 | 96.1                                       | 3.9             | 0.0              | 0.0               | 0.0          |

The ARCHITECT CA 19-9XR assay is used in conjunction with other clinical methods in the management of cancer patients.

It is recommended that each laboratory establish its own reference value for the population of interest.

\* Representative data; results in individual laboratories may vary from these data.

### Monitoring of Disease State in Patients Diagnosed with Pancreatic Cancer

Changes observed in serial CA 19-9 assay values when monitoring pancreatic cancer patients must be evaluated in conjunction with other clinical methods.

The effectiveness of the ARCHITECT CA 19-9XR assay as an aid in monitoring of disease state in pancreatic cancer patients was determined by assessing changes in levels of 1116-NS-19-9 reactive determinants in serial serum samples from 74 patients compared to changes in disease state. A study involving a total of 261 observations was performed with an average number of 3.5 observations per patient. In this study a significant change in levels of 1116-NS-19-9 reactive determinants was defined as at least a 14.0% increase in assay value (i.e., 2.5 times greater than the average of the assay's observed total %CV [5.6%]). A 14.0% change represents the minimum magnitude change between two serial ARCHITECT CA 19-9XR measurements that could not be attributed to assay variation or noise. Positive concordance between serial samples with at least a 14.0% increase in assay value and disease progression was found to be 48% (16/33). Negative concordance between serial samples with less than a 14.0% increase in assay value and no disease progression was found to be 64% (98/154). The overall concordance was found to be 61% (114/187). The following table presents the data in a 2 x 2 classification scheme\*.

| Change in Disease State per Sequential Pair               |             |                |       |
|---|-------------|----------------|-------|
| Change in the Level of 1116-NS-19-9 Reactive Determinants |             |                |       |
|   | Progression | No Progression | Total |
| ≥14.0%  | 16          | 56             | 72    |
| <14.0%  | 17          | 98             | 115   |
| Total   | 33          | 154            | 187   |



The following table provides the per patient distribution\*. Positive concordance between serial samples with at least a 14.0% increase in assay value and disease progression was found to be 68% (15/22). Negative concordance between serial samples with less than a 14.0% increase in assay value and no disease progression was found to be 69% (36/52). The overall concordance was found to be 69% (51/74).

| Change in Disease State per Patient                       |             |                |       |
|---|-------------|----------------|-------|
| Change in the Level of 1116-NS-19-9 Reactive Determinants | Progression | No Progression | Total |
| ≥14.0%  | 15          | 16             | 31    |
| <14.0%  | 7           | 36             | 43    |
| Total   | 22          | 52             | 74    |

\* Representative data; results in individual laboratories may vary from these data.

Below are examples of serial monitoring profiles for two patients with the disease state, ARCHITECT CA 19-9XR assay values, and the CA 19-9 RIA values.\* The disease states are:

- Progression from one collection to the next collection (Progression).
- No Change in disease state (Stable).
- Reduction in the signs and symptoms of the disease from one collection to the next (Responding).



\* Representative data; results in individual laboratories may vary from these data.

## SPECIFIC PERFORMANCE CHARACTERISTICS

### Precision

The ARCHITECT CA 19-9XR assay is designed to have an assay precision of ≤10% total CV.

A study was performed as described per the National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2.<sup>44</sup> Six samples were tested consisting of two panels of pooled serum (panels 1 and 2), one panel of serum to which 1116-NS-19-9 reactive determinants were added (panel 3), and the three ARCHITECT CA 19-9XR Controls. Testing was performed using two lots of reagents, in replicates of two at two separate times per day for 20

days on two separate instruments. Each reagent lot used a single calibration curve throughout the study. Data from this study are summarized below.\*

| Sample         | Reagent |            | n  | Mean Conc. (U/mL) | Within Run |     | Total |     |
|----------------|---------|------------|----|-------------------|------------|-----|-------|-----|
|                | Lot     | Instrument |    |                   | SD         | %CV | SD    | %CV |
| Panel 1        | 1       | 1          | 80 | 56.52             | 1.69       | 3.0 | 2.19  | 3.9 |
|                | 2       | 2          | 80 | 51.20             | 1.80       | 3.5 | 2.10  | 4.1 |
| Panel 2        | 1       | 1          | 80 | 311.49            | 7.22       | 2.3 | 10.72 | 3.4 |
|                | 2       | 2          | 80 | 288.82            | 9.14       | 3.2 | 11.23 | 3.9 |
| Panel 3        | 1       | 1          | 80 | 744.81            | 27.82      | 3.7 | 36.85 | 5.0 |
|                | 2       | 2          | 80 | 728.82            | 42.53      | 5.8 | 47.66 | 6.5 |
| Low Control    | 1       | 1          | 80 | 45.03             | 2.59       | 5.8 | 2.98  | 6.6 |
|                | 2       | 2          | 80 | 42.33             | 2.94       | 6.9 | 3.60  | 8.5 |
| Medium Control | 1       | 1          | 80 | 157.66            | 5.99       | 3.8 | 8.52  | 5.4 |
|                | 2       | 2          | 80 | 146.93            | 6.26       | 4.3 | 8.14  | 5.5 |
| High Control   | 1       | 1          | 80 | 781.68            | 44.76      | 5.7 | 49.87 | 6.4 |
|                | 2       | 2          | 80 | 781.42            | 62.10      | 8.0 | 65.28 | 8.4 |

\* Representative data; results in individual laboratories may vary from these data.

### Recovery

The ARCHITECT CA 19-9XR assay is designed to have a mean recovery of 100 ± 15% when 1116-NS-19-9 reactive determinants are added to serum samples.

A study was performed for the ARCHITECT CA 19-9XR assay based on guidance from Tietz Textbook of Clinical Chemistry.<sup>45</sup> Known concentrations of 1116-NS-19-9 reactive determinants were added to human serum samples. The concentration of 1116-NS-19-9 reactive determinants was determined using the ARCHITECT CA 19-9XR assay, and the resulting percent recovery was calculated. Representative data from this study are summarized in the table below.\*

| Sample | Endogenous Assay Value (U/mL) | 1116-NS-19-9 Reactive Determinants Added (U/mL) | Observed ARCHITECT CA 19-9XR Assay Value (U/mL) | % Recovery** |
|--------|-------------------------------|---|---|--------------|
| 1      | 46.50                         | 124.21  | 152.42  | 89           |
|        |                               | 629.91  | 645.00  | 95           |
| 2      | 28.96                         | 124.21  | 146.73  | 96           |
|        |                               | 629.91  | 598.93  | 91           |
| 3      | 38.42                         | 124.21  | 175.18  | 108          |
|        |                               | 629.91  | 652.12  | 98           |

Mean recovery across two separate spiked concentrations shown above = 96 %

$$** \% \text{ Recovery} = \frac{\text{Observed (U/mL)}}{\text{Endogenous Level (U/mL) + 1116-NS-19-9 Reactive Determinants Added (U/mL)}} \times 100$$

\* Representative data; results in individual laboratories may vary from these data.

### Dilution Linearity

The ARCHITECT CA 19-9XR assay is designed to have a mean recovery of 100 ± 15% of the expected result for diluted specimens.

A study was performed for the ARCHITECT CA 19-9XR assay modeled after the National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP6-P2.<sup>46</sup> Samples with known elevated 1116-NS-19-9 reactive determinant concentrations were diluted with ARCHITECT Multi-Assay Manual Diluent. The 1116-NS-19-9 reactive determinants concentration was determined for each dilution and the percent recovery was calculated. Representative data from this study are summarized below.\*

| Sample | Final Dilution Factor | Expected Value (U/mL) | Value Obtained (U/mL) | % Recovery** |
|--------|-----------------------|-----------------------|-----------------------|--------------|
| 1      | Undiluted             | 1024.55               | 1024.55               | -            |
|        | 1:2                   | 512.27                | 472.46                | 92           |
|        | 1:4                   | 256.14                | 264.26                | 103          |
|        | 1:5                   | 204.91                | 208.57                | 102          |
|        | 1:10                  | 102.45                | 108.94                | 106          |
|        | 1:20                  | 51.23                 | 54.33                 | 106          |
| 2      | Undiluted             | 1150.50               | 1150.50               | -            |
|        | 1:2                   | 575.25                | 551.62                | 96           |
|        | 1:4                   | 287.63                | 291.06                | 101          |
|        | 1:5                   | 230.10                | 253.65                | 110          |
|        | 1:10                  | 115.05                | 125.97                | 109          |
|        | 1:20                  | 57.53                 | 62.57                 | 109          |
| 3      | Undiluted             | 1028.25               | 1028.25               | -            |
|        | 1:2                   | 514.12                | 492.39                | 96           |
|        | 1:4                   | 257.06                | 290.24                | 113          |
|        | 1:5                   | 205.65                | 204.03                | 99           |
|        | 1:10                  | 102.82                | 120.76                | 117          |
|        | 1:20                  | 51.41                 | 57.25                 | 111          |

Mean recovery across the three diluted samples shown above = 105%

$$** \% \text{ Recovery} = \frac{\text{Values Obtained} \times \text{Dilution Factor}}{\text{Undiluted Expected Value}} \times 100$$

\* Representative data; results in individual laboratories may vary from these data.

#### Analytical Sensitivity

The analytical sensitivity of the ARCHITECT CA 19-9XR assay was calculated to be better than 2.00 U/mL (n = 18 runs, in replicates of 10).

Analytical sensitivity is defined as the concentration at two standard deviations from the ARCHITECT CA 19-9XR Calibrator A (0 U/mL), and represents the lowest measurable concentration of 1116-NS-19-9 reactive determinants that can be distinguished from zero.

#### Interference

The ARCHITECT CA 19-9XR assay is designed to have a mean recovery of 100 ± 12% in the presence of the chemotherapeutic agents listed below and elevated levels of bilirubin, hemoglobin, triglycerides, and total protein at the levels indicated.

A study based on guidance from the NCCLS Protocol EP7-A<sup>47</sup> was performed for the ARCHITECT CA 19-9XR assay. Specimens with 1116-NS-19-9 reactive determinant levels between 49.6 and 509.4 U/mL were supplemented with the following potentially interfering substances and chemotherapeutic agents.

#### POTENTIALLY INTERFERING SUBSTANCES

The average recovery observed during the study ranged from 91% to 102%.\*

| Substance     | Concentration |
|---------------|---------------|
| Bilirubin     | 22 mg/dL      |
| Hemoglobin    | 600 mg/dL     |
| Total Protein | 10 g/dL       |
| Triglycerides | 5100 mg/dL    |

\* Representative data; results in individual laboratories may vary from these data.

#### CHEMOTHERAPEUTIC AGENTS

The average recovery observed during the study ranged from 95% to 104%.\*

| Substance        | Concentration |
|------------------|---------------|
| 5-Fluorouracil   | 0.390 mg/mL   |
| Cisplatin        | 0.057 mg/mL   |
| Cyclophosphamide | 0.375 mg/mL   |
| Cytarabine       | 30 µg/mL      |
| Doxorubicin      | 40 µg/mL      |
| Gemcitabine      | 0.382 mg/mL   |
| Leucovorin       | 0.114 mg/mL   |
| Methotrexate     | 0.909 mg/mL   |
| Paclitaxel       | 0.067 mg/mL   |
| Streptozotocin   | 0.28 mg/mL    |
| Tamoxifen        | 2.28 µg/dL    |

\* Representative data; results in individual laboratories may vary from these data.

#### EVALUATION OF POTENTIALLY INTERFERING CLINICAL CONDITIONS

The ARCHITECT CA 19-9XR assay is designed to have a mean recovery of 100 ± 12% in the presence of HAMA and rheumatoid factor (RF).

The ARCHITECT CA 19-9XR assay was evaluated using specimens with HAMA and RF to further assess the clinical specificity. Five specimens positive for HAMA and five specimens positive for RF were evaluated for % recovery with 1116-NS-19-9 reactive determinants spiked into each specimen at 35 and 250 U/mL. Mean percent recovery results are summarized in the following table.\*

| Clinical Condition | Number of Specimens | Mean % Recovery |
|--------------------|---------------------|-----------------|
| HAMA               | 10                  | 93              |
| RF                 | 10                  | 93              |

\* Representative data; results in individual laboratories may vary from these data.

#### Carryover

No significant carryover (less than 2.00 U/mL in CA19-9XR Calibrator A\*) was observed for the ARCHITECT CA 19-9XR assay when a sample containing up to 320,000 U/mL of 1116-NS-19-9 reactive determinants was assayed.

\* Representative data; results in individual laboratories may vary from these data.

#### High Dose Hook

No high dose hook effect was observed for the ARCHITECT CA 19-9XR assay when samples containing up to 1,750,000 U/mL\* of 1116-NS-19-9 reactive determinants were assayed. High dose hook is a phenomenon whereby very high level specimens may falsely read within the dynamic range of the assay.






\* Representative data; results in individual laboratories may vary from these data.

## BIBLIOGRAPHY

- Hotakainen K, Tanner P, Alftan H, et al. Comparison of three immunoassays for CA 19-9. *Clin Chim Acta*. 2009;400:123-127.
- Passerinin R, Cassatella M, Boveri S, et al. Routine Testing and Comparison of Two Automated Immunoassays in a Reference Oncology Center. *Am J Clin Pathol*. 2012;138:281-287.
- Uhlenbruck G, Höller U, Heising J, et al. Sialylated Le<sup>a</sup> Blood Group Substances Detected by the Monoclonal Antibody CA 19-9 in Human Seminal Plasma and Other Organs. *Urol Res* (Germany, West) 1985;13(15):223-226.
- Herlyn M, Stepkowski Z, Herlyn D, et al. Colorectal Carcinoma-Specific Antigen: Detection by Means of Monoclonal Antibodies. *Proc Natl Acad Sci USA* 1979;76:1438-1442.
- Magnani JL, Brockhaus M, Smith DF, et al. A Monosialoganglioside is a Monoclonal Antibody-Defined Antigen of Colon Carcinoma. *Science* 1981;212:55-56.
- Magnani J, Nilsson B, Brockhaus M, et al. The Antigen of a Tumor-Specific Monoclonal Antibody is a Ganglioside Containing Sialylated Lacto-N-Fucopentaose II. *Fed Proc* 1982;41:898.
- Magnani JL, Stepkowski Z, Koprowski H, et al. Identification of the Gastrointestinal and Pancreatic Cancer-Associated Antigen Detected by Monoclonal Antibody 19-9 in the Sera of Patients as a Mucin. *Cancer Res* 1983;43:5489-5492.
- Del Villano BC, Brennan S, Brock P, et al. Radioimmunoassay for a Monoclonal Antibody-Defined Tumor Marker, CA 19-9. *Clin Chem* 1983;29:549-552.
- Steinberg WM, Gelfand R, Anderson KK, et al. Comparison of the Sensitivity and Specificity of the CA 19-9 and Carcinoembryonic Antigen Assays in Detecting Cancer of the Pancreas. *Gastroent* 1986;90:343-349.
- Ritts RE Jr, Del Villano BC, Go VLW, et al. Initial Clinical Evaluation of an Immunoradiometric Assay for CA 19-9 Using the NCI Serum Bank. *Int J Cancer* 1984;33:339-345.
- Jalanko H, Kuusela P, Roberts P, et al. Comparison of a New Tumor Marker, CA 19-9, with Alpha-Fetoprotein and Carcinoembryonic Antigen in Patients with Upper Gastrointestinal Diseases. *J Clin Pathol* 1984;37:218-222.
- Gupta MK, Arciaga R, Bocci L, et al. Measurement of a Monoclonal Antibody-Defined Antigen (CA 19-9) in the Sera of Patients with Malignant and Nonmalignant Diseases. Comparison with Carcinoembryonic Antigen. *Cancer* 1985;56(2):277-283.
- Andriulli A, Gindro T, Piantino P, et al. Prospective Evaluation of the Diagnostic Efficacy of CA 19-9 Assay as a Marker for Gastrointestinal Cancers. *Digestion* 1986;33(1):26-33.
- Patai Á, Héber S, Döbrönte Z, et al. Diagnostic Values of CA 19-9 and CEA in Gastrointestinal Diseases. *Orv Hetil* 1992;133(21):1301-1304,1307.
- Ritts RE and Pitt HA. CA 19-9 in Pancreatic Cancer. *Surgical Oncology Clinics of North America* 1998;7:93-101.
- Duffy MJ, van Dalen A, Haglund C, et al. Clinical Utility of Biochemical Markers in Colorectal Cancer: European Group on Tumor Markers (EGTM) Guidelines. *Eur J Cancer* 2003;39:718-727.
- Cerwenka H, Aigner R, Quehenberger F, et al. Preoperative Differential Diagnosis of Benign and Malignant Pancreatic Lesions - The Value of Pancreatic Secretory Trypsin Inhibitor, Procarboxypeptidase B, CA 19-9 and CEA. *Hepato-Gastroenterology* 1997;44(16):1117-1121.
- Von Ritter C, Eder MI, Stieber P, et al. Biliary Mucin Secreted by Cultured Human Gallbladder Epithelial Cells Carries the Epitope of CA 19-9. *Anticancer Res* 1997;17(4B):2931-2934.
- Adachi Y, Iso Y, Moriyama M, et al. Increased Serum CA 19-9 in Patients with Xanthogranulomatous Cholecystitis. *Hepta-Gastroenterology* 1998;45:77-80.
- Maestranzi S, Przemioslo R, Mitchell H, et al. The Effect of Benign and Malignant Liver Disease on the Tumor Markers CA 19-9 and CEA. *Ann Clin Biochem* 1998;35:99-103.
- Duffy MJ, O'Sullivan F, McDonnell TJ, et al. Increased Concentrations of the Antigen CA 19-9 in Serum of Cystic Fibrosis Patients. *Clin Chem* 1985;31:1245-1246 (Letter).
- Kane RE, Penny J, Walker K, et al. Changes in the CA 19-9 Antigen and Lewis Blood Group with Pulmonary Disease Severity in Cystic Fibrosis. *Pediatr Pulmonol* 1992;12(4):221-226.
- Wu JT and Chang J. Chromatographic Characterization of CA 19-9 Molecules from Cystic Fibrosis and Pancreatic Carcinoma. *J Clin Lab Anal* 1992;6(4):209-215.
- Wu JT, Olsen J, Walker K, Tumor Markers CA 19-9 and CA 195 Are Also Useful as Markers for Cystic Fibrosis. *J Clin Lab Anal* 1992;5(3):151-161.
- Staab HJ, Brümmendorf T, Hornung A, et al. The Clinical Validity of Circulating Tumor-Associated Antigens CEA and CA 19-9 in Primary Diagnosis and Follow-up of Patients with Gastrointestinal Malignancies. *Klin Wochenschr* 1985;63(3):106-115.
- Ychou M, Tuszynski T, Pignon J-P, et al. Gastric Carcinoma: Comparison between CEA and CA 19-9 for diagnosis of recurrence after gastrectomy. *Gastroenterol Clin Biol* 1992;16(11):848-852.
- Grem J. The Prognostic Importance of Tumor Markers in Adenocarcinomas of the Gastrointestinal Tract. *Curr Opin Oncol* 1997;9(4):380-387.
- Gärtner U, Scheulen ME, Conradt C, et al. Value of Tumor-Associated Antigen CA 72-4 Compared with CEA and CA 19-9 in the Followup of Patients Operated for Gastric Carcinoma. *Dtsch med Wschr* 1998;123:69-73.
- Willet CG, Daly WJ, Warshaw AL, CA 19-9 is an Index of Response to Neoadjuvant Chemoradiation Therapy in Pancreatic Cancer. *Am J Surg* 1996;172(4):350-352.
- Filella X, Molina R, Grau JJ, et al. Prognostic Value of CA 19.9 Levels in Colorectal Cancer. *Ann Surg* 1992;216(1):55-59.
- Kouri M, Pyrhönen S, Kuusela P, Elevated CA 19-9 as the Most Significant Prognostic Factor in Advanced Colorectal Carcinoma. *J Surg Oncol* 1992;49(2):78-85.
- Gebauer G, Müller-Ruchholtz W. Tumor Marker Concentrations in Normal and Malignant Tissues of Colorectal Cancer Patients and Their Prognostic Relevance. *Anticancer Res* 1997;17(4a):2731-2734.
- Reiter W, Stieber P, Reuter C, et al. Preoperative Serum Levels of CEA and CA 19-9 and Their Prognostic Significance in Colorectal Carcinoma. *Anticancer Res* 1997;17(4B):2935-2938.
- Gogas H, Lofts FJ, Evans TRJ, et al. Are Serial Measurements of CA 19-9 Useful in Predicting Response to Chemotherapy in Patients with Inoperable Adenocarcinoma of the Pancreas? *Br J Cancer* 1998;77:325-328.
- Safi F, Schlosser W, Falkenreck S, et al. Prognostic Value of CA 19-9 Serum Course in Pancreatic Cancer. *Hepato-Gastroenterol* 1998;45:253-259.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
- Tempero MA, Uchida E, Takasaki H, et al. Relationship of Carbohydrate Antigen 19-9 and Lewis Antigens in Pancreatic Cancer. *Cancer Res* 1987;47:5501-5503.
- National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
- Tietz NW. *Tietz Textbook of Clinical Chemistry*. 3rd Edition. Philadelphia, PA: B. Saunders Company, 1999;330.
- National Committee for Clinical Laboratory Standards. *Evaluation of the Linearity of Quantitative Analytical Methods; Proposed Guideline—Second Edition*. NCCLS Document EP6-P2. Wayne, PA: NCCLS, 2001.
- National Committee for Clinical Laboratory Standards (NCCLS). *Interference Testing in Clinical Chemistry; Approved Guideline*. NCCLS Document EP7-A. Wayne, PA: NCCLS; 2002.



## ■ Key to Symbols

|   |   |
|---|---|
|  | Consult instructions for use  |
|  | Manufacturer  |
|  | Sufficient for  |
|  | Temperature limitation  |
|  | Use by/Expiration date  |
| <b>CONJUGATE</b>  | Conjugate   |
| <b>CONTAINS: AZIDE</b>  | Contains Sodium Azide. Contact with acids liberates very toxic gas. |
| <b>CONTROL NO.</b>  | Control Number  |
| <b>ECO HAZARD</b>   | Ecological hazard   |
| <b>IVD</b>  | <i>In Vitro</i> Diagnostic Medical Device                           |
| <b>LOT</b>  | Lot Number  |
| <b>MICROPARTICLES</b>   | Microparticles  |
| <b>MULTI-ASSAY MANUAL DILUENT</b>   | Multi-Assay Manual Diluent  |
| <b>PRE-TRIGGER SOLUTION</b>   | Pre-Trigger Solution  |
| <b>PRODUCED FOR ABBOTT BY</b>   | Produced for Abbott by  |
| <b>PRODUCT OF USA</b>   | Product of USA  |
| <b>REACTION VESSELS</b>   | Reaction Vessels  |
| <b>REAGENT LOT</b>  | Reagent Lot   |
| <b>REF</b>  | List Number   |
| <b>REPLACEMENT CAPS</b>   | Replacement Caps  |
| <b>SAMPLE CUPS</b>  | Sample Cups   |
| <b>SEPTUM</b>   | Septum  |
| <b>SN</b>   | Serial number   |
| <b>TRIGGER SOLUTION</b>   | Trigger Solution  |
| <b>WARNING: SENSITIZER</b>  | Warning: May cause an allergic reaction.                            |
| <b>WASH BUFFER</b>  | Wash Buffer   |

ARCHITECT and Chemiflex are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.



Abbott GmbH & Co. KG  
Max-Planck-Ring 2  
65205 Wiesbaden  
Germany  
+49-6122-580



### **PRODUCED FOR ABBOTT BY**

Fujirebio Diagnostics Inc., Malvern, PA 19355 USA

**Customer Service: Contact your local representative or find country-specific contact information on [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)**

Revised November 2016.

©2005, 2016 Abbott Laboratories

