



**ARCHITECT**  
**Sirolimus Whole Blood**  
**Precipitation Reagent**

REF 1L76-55



093374.00 REV 000 / R07  
H1L760



**en**  
**Sirolimus**  
**1L76**  
**093374.00 REV 000 / R07**  
**H1L760**

Read Highlighted Changes: Revised June 2017.

**INTENDED USE**

The ARCHITECT Sirolimus Whole Blood Precipitation Reagent is for the extraction of sirolimus from samples (human whole blood patient specimens, controls and ARCHITECT Sirolimus Calibrators) to be tested on the ARCHITECT iSystem.

Refer to the ARCHITECT Sirolimus reagent package insert for pretreatment instructions.

**CONTENTS**

1 Bottle (30 mL) ARCHITECT Sirolimus Whole Blood Precipitation Reagent containing zinc sulfate solution in DMSO and ethylene glycol.

**PRECAUTIONS**

- **IVD**
- For *In Vitro* Diagnostic Use
- **CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.<sup>1-4</sup>


The following warnings and precautions apply to: <b>WHOLE BLOOD PRECIPITATION REAGENT</b>	
<b>WARNING</b>	Contains ethylene glycol, dimethyl sulfoxide and zinc sulfate.
H373	May cause damage to organs through prolonged or repeated exposure.
H319	Causes serious eye irritation.
H412	Harmful to aquatic life with long lasting effects.
<b>Prevention</b>	
P260	Do not breathe mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.
P273	Avoid release to the environment.
<b>Response</b>	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice / attention.
P314	Get medical advice / attention if you feel unwell.
<b>Disposal</b>	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at [www.abbottdiagnostics.com](http://www.abbottdiagnostics.com) or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

**STORAGE**

- Reagents are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.



**BIBLIOGRAPHY**

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

**Key to Symbols**

	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
<b>EC REP</b>	Authorized Representative in the European Community
<b>INFORMATION FOR USA ONLY</b>	Information needed for United States of America only
<b>IVD</b>	<i>In Vitro</i> Diagnostic Medical Device
<b>LOT</b>	Lot Number
<b>PRODUCED FOR ABBOTT BY</b>	Produced for Abbott by
<b>PRODUCT OF USA</b>	Product of USA
<b>REF</b>	List Number
<b>WHOLE BLOOD PRECIPITATION REAGENT</b>	Whole Blood Precipitation Reagent

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.



Abbott Laboratories  
Diagnostics Division  
Abbott Park, IL 60064  
USA



Abbott GmbH & Co. KG  
Max-Planck-Ring 2  
65205 Wiesbaden  
Germany  
+49-6122-580



**PRODUCED FOR ABBOTT BY**

Fujirebio Diagnostics Inc., Seguin, TX 78155 USA

**Customer Service: Contact your local representative  
or find country-specific contact information on  
[www.abbottdiagnostics.com](http://www.abbottdiagnostics.com)**

Revised June 2017.

©2006, 2017 Abbott Laboratories

