



HI-TORQUE BALANCE MIDDLEWEIGHT™ GUIDE WIRES

TRUSTED DESIGN

NITINOL DISTAL CORE

for durability

SHAPING RIBBON

for a soft atraumatic tip

HYDROPHILIC OR HYDROPHOBIC COATING

for choice of lubricity level



Data on file at Abbott.

**Information contained herein for DISTRIBUTION outside the U.S. ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**


©2021 Abbott. All rights reserved. MAT-2102011 v1.0

HI-TORQUE BALANCE MIDDLEWEIGHT™ GUIDE WIRE ORDERING INFORMATION



HI-TORQUE BALANCE MIDDLEWEIGHT™

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)	DISTAL COATING
HI-TORQUE BALANCE MIDDLEWEIGHT™	1001780-HC	190	Straight	3	0.6	Hydrophilic
HI-TORQUE BALANCE MIDDLEWEIGHT™	1001780J-HC	190	J	3	0.6	Hydrophilic
HI-TORQUE BALANCE MIDDLEWEIGHT™	1001782-HC	300	Straight	3	0.6	Hydrophilic
HI-TORQUE BALANCE MIDDLEWEIGHT™	1001782J-HC	300	J	3	0.6	Hydrophilic
HI-TORQUE BALANCE MIDDLEWEIGHT™	1001780	190	Straight	3	0.6	Hydrophobic
HI-TORQUE BALANCE MIDDLEWEIGHT™	1001780J	190	J	3	0.6	Hydrophobic
HI-TORQUE BALANCE MIDDLEWEIGHT™	1001782	300	Straight	3	0.6	Hydrophobic
HI-TORQUE BALANCE MIDDLEWEIGHT™	1001782J	300	J	3	0.6	Hydrophobic



*Tests performed by and data on file at Abbott.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

Information contained herein for DISTRIBUTION outside the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Data on file at Abbott.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

©2021 Abbott. All rights reserved. MAT-2102011 v1.0