



HI-TORQUE BALANCE HEAVYWEIGHT™

Guide Wires

TRUSTED DURABILITY WITH ADDED SUPPORT

NITINOL DISTAL CORE

for durability

ADDED DISTAL SUPPORT

for device delivery and vessel straightening

SHAPING RIBBON

for a soft atraumatic tip



Data on file at Abbott.

**Information contained herein for DISTRIBUTION outside the U.S. ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**


©2021 Abbott. All rights reserved. MAT-2102010 v1.0

HI-TORQUE BALANCE HEAVYWEIGHT™ GUIDE WIRE ORDERING INFORMATION



HI-TORQUE BALANCE HEAVYWEIGHT™

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE BALANCE HEAVYWEIGHT™	1000462H	190	Straight	4.5	0.9
HI-TORQUE BALANCE HEAVYWEIGHT™	1000462HJ	190	J	4.5	0.9
HI-TORQUE BALANCE HEAVYWEIGHT™	1000463H	300	Straight	4.5	0.9
HI-TORQUE BALANCE HEAVYWEIGHT™	1000463HJ	300	J	4.5	0.9



*Tests performed by and data on file at Abbott.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

Information contained herein for DISTRIBUTION outside the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Data on file at Abbott.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

©2021 Abbott. All rights reserved. MAT-2102010 v1.0

