

# HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL™ & HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™ GUIDE WIRES

## TRUSTED DESIGN

### **NITINOL DISTAL CORE**

for durability

### **SHAPING RIBBON**

for a soft atraumatic tip

### **INTERMEDIATE POLYMER SLEEVE**

for smooth tracking



Data on file at Abbott.

**Information contained herein for DISTRIBUTION outside the U.S. ONLY.  
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

©2021 Abbott. All rights reserved. MAT-2102012 v1.0

# HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL™ & HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™ GUIDE WIRE ORDERING INFORMATION



## HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL™

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL™	1009660	190	Straight	3	0.5
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL™	1009660J	190	J	3	0.5
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL™	1009661	300	Straight	3	0.5
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL™	1009661J	300	J	3	0.5

## HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)	
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™	1009664	190	Straight	3	0.6	WITH MARKERS
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™	1009664J	190	J	3	0.6	
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™	1009665	300	Straight	3	0.6	WITHOUT MARKERS
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™	1009665J	300	J	3	0.6	
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™	1009666	190	Straight	3	0.6	WITHOUT MARKERS
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™	1009666J	190	J	3	0.6	
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™	1009667	300	Straight	3	0.6	
HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™	1009667J	300	J	3	0.6	



\*Tests performed by and data on file at Abbott.

^Optional for HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™ Guide Wire.

\*\*HI-TORQUE BALANCE MIDDLEWEIGHT UNIVERSAL II™ Guide Wire has an additional silicone-based hydrophobic coating.

**CAUTION:** This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at [vascular.eifu.abbott](http://vascular.eifu.abbott) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

**Information contained herein for DISTRIBUTION outside the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Data on file at Abbott.

### Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

[www.cardiovascular.abbott](http://www.cardiovascular.abbott)

©2021 Abbott. All rights reserved. MAT-2102012 v1.0