

# HI-TORQUE POWERTURN™

Guide Wires

## THE POWER OF TIP CONTROL

### **CORE-TO-TIP DESIGN**


for steering and tip control

### **PARABOLIC CORE GRIND**

for vessel tracking and steering

### **3 SUPPORT LEVEL OFFERINGS**

to meet a range of clinical needs in workhorse cases



Data on file at Abbott.

**Information contained herein for DISTRIBUTION outside the U.S. ONLY.  
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

©2021 Abbott. All rights reserved. MAT-2102004 v1.0

# HI-TORQUE POWERTURN™ GUIDE WIRE ORDERING INFORMATION



## HI-TORQUE POWERTURN™ ULTRAFLEX

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE POWERTURN™ ULTRAFLEX	1044590	190	Straight	3	1.0
HI-TORQUE POWERTURN™ ULTRAFLEX	1044590J	190	J	3	1.0
HI-TORQUE POWERTURN™ ULTRAFLEX	1044591	300	Straight	3	1.0
HI-TORQUE POWERTURN™ ULTRAFLEX	1044591J	300	J	3	1.0

## HI-TORQUE POWERTURN™ FLEX

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE POWERTURN™ FLEX	1044592	190	Straight	3	1.0
HI-TORQUE POWERTURN™ FLEX	1044592J	190	J	3	1.0
HI-TORQUE POWERTURN™ FLEX	1044593	300	Straight	3	1.0
HI-TORQUE POWERTURN™ FLEX	1044593J	300	J	3	1.0

## HI-TORQUE POWERTURN™

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE POWERTURN™	1044594	190	Straight	3	0.9
HI-TORQUE POWERTURN™	1044594J	190	J	3	0.9
HI-TORQUE POWERTURN™	1044595	300	Straight	3	0.9
HI-TORQUE POWERTURN™	1044595J	300	J	3	0.9

\*Tests performed by and data on file at Abbott.

**CAUTION:** This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at [vascular.eifu.abbott](http://vascular.eifu.abbott) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

**Information contained herein for DISTRIBUTION outside the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Data on file at Abbott.

### Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

[www.cardiovascular.abbott](http://www.cardiovascular.abbott)

©2021 Abbott. All rights reserved. MAT-2102004 v1.0

