

HI-TORQUE TURNTRAC™

Guide Wires

1:1 TORQUE RESPONSE¹ FOR CHALLENGING WORKHORSE CASES

NITINOL TIP

for durability

CORE-TO-TIP DESIGN

for excellent tip control and steering*

PARABOLIC CORE GRIND

for vessel tracking and torque response

TIP COATING OPTIONS


Uncoated tip for increased tactile feel

Optional coated tip for lubricity

1. Tóth GG, Yamane M, Heyndrickx GR. Heart 2015;101:645–652.

*Tests performed by and data on file at Abbott.

**Information contained herein for DISTRIBUTION outside the U.S. ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**



HI-TORQUE TURNTRAC™ GUIDE WIRE ORDERING INFORMATION



HI-TORQUE TURNTRAC™ FLEX


PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	TIP COATING	TIP LOAD*(g)
HI-TORQUE TURNTRAC™ FLEX	1020009	190	Straight	Coated	0.6
HI-TORQUE TURNTRAC™ FLEX	1020010	190	Straight	Uncoated	0.6
HI-TORQUE TURNTRAC™ FLEX	1020011	300	Straight	Coated	0.6
HI-TORQUE TURNTRAC™ FLEX	1020012	300	Straight	Uncoated	0.6
HI-TORQUE TURNTRAC™ FLEX	1020009J	190	J	Coated	0.6
HI-TORQUE TURNTRAC™ FLEX	1020010J	190	J	Uncoated	0.6

HI-TORQUE TURNTRAC™

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	TIP COATING	TIP LOAD*(g)
HI-TORQUE TURNTRAC™	1020013	190	Straight	Coated	0.8
HI-TORQUE TURNTRAC™	1020014	190	Straight	Uncoated	0.8
HI-TORQUE TURNTRAC™	1020015	300	Straight	Coated	0.8
HI-TORQUE TURNTRAC™	1020016	300	Straight	Uncoated	0.8
HI-TORQUE TURNTRAC™	1020013J	190	J	Coated	0.8
HI-TORQUE TURNTRAC™	1020014J	190	J	Uncoated	0.8



**TIP-SHAPING TOOL
INCLUDED WITH EVERY WIRE**
for convenient tip-shaping



*Tests performed by and data on file at Abbott.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

Information contained herein for DISTRIBUTION outside the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Data on file at Abbott.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

©2021 Abbott. All rights reserved. MAT-2102013 v1.0

