

HI-TORQUE WHISPER™

Guide Wires

DESIGNED FOR SMOOTH
DEVICE DELIVERY IN HIGHLY
TORTUOUS ANATOMY

CORE-TO-TIP DESIGN

for excellent torque response¹

FULL POLYMER COVER

for smooth vessel tracking

POLYMER COVER AND HYDROPHILIC COATING

for outstanding lubricity¹

3 SUPPORT LEVEL OPTIONS

for crossing and device delivery

Data on file at Abbott.

1. Tóth GG, Yamane M, Heyndrickx GR. Heart 2015;101:645–652.

Information contained herein for DISTRIBUTION outside the U.S. ONLY.

Check the regulatory status of the device in areas where CE marking is not the regulation in force.

©2021 Abbott. All rights reserved. MAT-2102008 v1.0



HI-TORQUE WHISPER™ GUIDE WIRE ORDERING INFORMATION



HI-TORQUE WHISPER™ ES


PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE WHISPER™ ES	1011834H	190	Straight	3	1.3
HI-TORQUE WHISPER™ ES	1011834HJ	190	J	3	1.3
HI-TORQUE WHISPER™ ES	1011835H	300	Straight	3	1.3
HI-TORQUE WHISPER™ ES	1011835HJ	300	J	3	1.3

HI-TORQUE WHISPER™ MS

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE WHISPER™ MS	1005357H	190	Straight	3	1.0
HI-TORQUE WHISPER™ MS	1005357HJ	190	J	3	1.0
HI-TORQUE WHISPER™ MS	1005359H	300	Straight	3	1.0
HI-TORQUE WHISPER™ MS	1005359HJ	300	J	3	1.0

HI-TORQUE WHISPER™ LS

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE WHISPER™ LS	1005351H	190	Straight	3	1.0
HI-TORQUE WHISPER™ LS	1005351HJ	190	J	3	1.0
HI-TORQUE WHISPER™ LS	1005353H	300	Straight	3	1.0
HI-TORQUE WHISPER™ LS	1005353HJ	300	J	3	1.0



*Tests performed by and data on file at Abbott.

1. Tóth GG, Yamane M, Heyndrickx GR. Heart 2015;101:645–652.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

Information contained herein for DISTRIBUTION outside the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Data on file at Abbott.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

©2021 Abbott. All rights reserved. MAT-2102008 v1.0