

CORONARY

LEADING THE WAY IN GUIDE WIRE TECHNOLOGY

Communication **HI-TORQUE INFILTRAC[™]**

Guide Wire Family

A CTO INDICATED PENETRATION GUIDE WIRE

Information contained herein for DISTRIBUTION outside of the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

HI-TORQUE INFILTRAC[™] GUIDE WIRE'S STATE-OF-THE-ART TIP TECHNOLOGY^{*} IS DESIGNED FOR PENETRATING RESISTANT CTO CAPS^{*}

Unique Micro-Textured Tip for Cap Engagement[^]



- **HI-TORQUE INFILTRAC**[™] Guide Wire's unique micro-textured tip provides excellent traction[^]
- **HI-TORQUE INFILTRAC**[™] Guide Wire's uncoated tip delivers exceptional tactile feel*

Tapered Tip for High Penetration Powers[^]

TWO PENETRATION POWER OPTIONS FOR ANTEGRADE WIRE ESCALATION

- HI-TORQUE INFILTRAC[™] Guide Wire = 167 kg/in² (11 g tip load)
- HI-TORQUE INFILTRAC[™] PLUS Guide Wire = 224 kg/in² (14 g tip load)

.009 in. Tip Diameter

*Data on file at Abbott. ^Tests performed by and data on file at Abbott.

Information contained herein for DISTRIBUTION outside of the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.


Pre-shaped Micro-J Tip for Enhanced Steering*

1 mm from tip



• HI-TORQUE INFILTRAC[™] Guide Wire's pre-shaped micro-j tip provides steering capabilities while still maintaining high penetration powers^{*}

HI-TORQUE INFILTRAC[™] Guide Wire Combines Unique Tip Technology with Proven Guide Wire Design Features^{*}



SPECIFICATION	HI-TORQUE INFILTRAC™ GUIDE WIRE FAMILY			
Core Material	DURASTEEL Stainless Steel			
Cover	None (bare coils)			
Distal Coating Type	Hydrophilic			
Tip Coating	Uncoated			
Tip Shape	Micro-J			
Tip Diameter (in)	0.009 (tapered)			
Radiopaque Coil Length (cm)	3			
Wire Length (cm)	190, 300			

*Data on file at Abbott.

'Tests performed by and data on file at Abbott.

Information contained herein for DISTRIBUTION outside of the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Ordering Information*

DESCRIPTION	NON-CE" PART NUMBER	CE PART NUMBER	TIP LOAD^ (g)	PENETRATION POWER^ (kg/in ²)	WIRE LENGTH (cm)	TIP DIAMETER (in)	TIP SHAPE
HI-TORQUE INFILTRAC [™] Guide Wire	1030001J	1030005J	11	167	190	0.009	Micro-J
	1030002J	1030006J	11	167	300	0.009	Micro-J
HI-TORQUE INFILTRAC [™] PLUS Guide Wire	1030003J	1030007J	14	224	190	0.009	Micro-J
	1030004J	1030008J	14	224	300	0.009	Micro-J

*Data on file at Abbott.

**Same as U.S. FDA part numbers applicable to select OUS countries. ^Tests performed by and data on file at Abbott.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at *vascular.eifu.abbott* or at *medical.abbott/manuals* for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for DISTRIBUTION outside of the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

Abbott International BVBA Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11 ™ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

@2021 Abbott. All rights reserved. MAT-2100873 v1.0

