

CAROTID

EMBOSHIELD NAV6TM

Embolic Protection System

NAVIGATION. ACCESS. VISIBILITY.

Navigate with a new level of control, confidence and accuracy

Navigation


• Distal access, workhorse and support BareWires designed specifically for carotid anatomies to ensure navigational success through tortuous anatomy and challenging arches

Access

 BareWire, fixed and telescope filter delivery options provide a customized approach to access and support diverse vascular anatomies

Visibility

 Radiopaque circumferential frame conveys vessel apposition and deployment status



BAREWIRE FILTER DELIVERY WIRES

Distal access

A soft barewire with extra flexibility and torqueability for successful access through tortuous anatomy.



Workhorse

Provides optimal flexibility and suport for superb access through most anatomies.




Support


Designed to provide additional support for delivery of filter in challenging arches.



FILTER VISIBILITY



Fully deployed



Fully retracted

ORDERING INFORMATION

EMBOSHIELD NAV ^{6TM} EMBOLIC PROTECTION SYSTEM*						
ORDER NUMBER	FILTER DIAMETER EXPANDS (MM)	RECOMMENDED VESSEL DIAMETER (MM)		MIN. SIZE OF SHEATH/ GUIDE CATHETER	MIN. INTERNAL DIAMETER SHEATH/GUIDE CATHETER (IN)	
22442-19	5	2.5 – 4.8	190	5/6	0.070	
22443-19	7.2	4-7	190	5/6	0.070	

^{*}EMBOSHIELD NAV^{6™} consists of 2 packages:

Package	Contents
Filter system	(1) RX application catheter, (1) filter element, (1) BareWire filter delivery wire
	(workhorse, 190 cm), (1) loading funnel, (1) introducer, (1) torquer, (1) syringe
RX retrieval catheter	with rinse tip (1) RX retrieval catheter, (1) torquer, (1) syringe with rinse tip

BAREWIRE DELIVERY WIRES (0.014")**					
ORDER NUMBER	WIRE TYPE	WIRE LENGTH (CM)			
22444-19	Distal Access	190			
22445-19 22445-31	Workhorse	190 315			
22446-19	Support	190			

 $^{**\}mbox{BareWires}$ are available in packs of 5

Caution: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use provided inside the product carton (when available), at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. **Information contained herein for DISTRIBUTION in Europe, Middle East and Africa ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Tests performed by and data on file at Abbott. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

Abbott Vascular International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

 $^{\mbox{\tiny TM}}$ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

©2019 Abbott. All rights reserved. 2-EH-1-9299-01 03-2019

