


Don't Compromise. Take Command.

Hybrid workhorse wires designed for below-the-knee success.



Excellent torque and support


- Stainless steel body and a nitinol tip that provides the optimal balance of control and durability
- Transitionless fusion weld for smooth device compatibility
- Hydrophilic coating on a radiopaque polymer for excellent lubricity and visibility
- Core-to-tip design enables 1:1 torque over the length of the wire*

Superb flexibility & tip durability

- Shapeable nitinol tip that retains shape throughout an entire procedure
- 3 cm distal platinum radiopaque coils for excellent visibility*

Wire and tip options

- Tip load options of 2.8 g and 3.5 g
- Wire lengths of 190 cm, 250 cm and 300 cm



Tests performed by and data on file at Abbott.

*Tóth GG, Yamane M, Heyndrickx GR. Heart 2015;101:645-652

Information contained herein for DISTRIBUTION outside the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

©2021 Abbott. All rights reserved. MAT-2010684 v1.0

Hi-Torque Command™

.014 Workhorse Guide Wires

Ordering Information

DESCRIPTION	WIRE LENGTH	PART NUMBER	TIP STYLE	UNITS PER PACKAGE
Hi-Torque Command™	190 cm	2078172	Core-To-Tip	5
Hi-Torque Command™	250 cm	2078176	Core-To-Tip	5
Hi-Torque Command™	300 cm	2078173	Core-To-Tip	5
Hi-Torque Command™ ES	190 cm	2078174	Core-To-Tip	5
Hi-Torque Command™ ES	250 cm	2078177	Core-To-Tip	5
Hi-Torque Command™ ES	300 cm	2078175	Core-To-Tip	5

Data on file at Abbott

CAUTION: These products are intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at www.vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

Information contained herein for DISTRIBUTION outside the U.S. only. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

©2021 Abbott. All rights reserved. MAT-2010684 v1.0

