

Hi-Torque Winn™ Family

.014 Specialty Guide Wires

PERIPHERAL

.014 | GUIDE WIRE

Feel. Cross. Winn.

A family of specialty guide wires for crossing challenging lesions.



Feel

- The excellent tactile feedback and deliverability¹ are due to the combination of:
 - Intermediate polymer sleeve
 - Exposed shapeable tip coils (distal 5 mm)

Cross

- Core-to-tip wire design, with transitionless durasteel core taper, enables pushability as well as outstanding torque over the length of the wire¹
- Hydrophilic coating provides excellent lubricity and ease of wire advancement¹

Winn

- Three distinct guide wires that offer a selection of tip loads and penetration power
- Tip load options of 5.1 g, 11.3 g, and 14.3 g
- Wire lengths of 190 cm and 300 cm

Ordering Information

PRODUCT	WIRE LENGTH	STOCK NUMBER	TIP STYLE	TIP LOAD (grams)
Hi-Torque Winn™ 40	190 cm	1012466	Core-To-Tip	5.1
Hi-Torque Winn™ 40	300 cm	1012467	Core-To-Tip	5.1
Hi-Torque Winn™ 80	190 cm	1012468	Core-To-Tip	11.3
Hi-Torque Winn™ 80	300 cm	1012469	Core-To-Tip	11.3
Hi-Torque Winn™ 200T	190 cm	1012474	Core-To-Tip	14.3
Hi-Torque Winn™ 200T	300 cm	1012475	Core-To-Tip	14.3

¹ Tóth GG, Yamane M, Heyndrickx GR. Heart 2015;101:645–652. Test(s) performed by and data on file at Abbott

Information contained herein for DISTRIBUTION outside of the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

©2021 Abbott. All rights reserved. MAT-2103475 v1.0



CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at *vascular.eifu.abbott* or at *medical.abbott/manuals* for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Information contained herein for DISTRIBUTION outside of the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

TM Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

©2021 Abbott. All rights reserved. MAT-2103475 v1.0

