

TREK™

Coronary Balloon Dilatation Catheter

Designed for challenging environments

Available in RX and OTW

CHALLENGING ANATOMY

Transitionless Tip

For smooth, gradual transition, providing ultimate control for lesion crossing and stent strut deflection



SUPERB TRACKABILITY

Multi-Layer CrossFlex²™ Balloon* and Dual Tungsten Markers

Enables smooth, precise tracking while providing for controlled and confident placement



TREK™ OTW

Tapered Proximal Shaft† With Low Profile Distal Shaft††

Provides a balance of pushability and flexibility

OUTSTANDING DELIVERABILITY

Slim Seal™ Technology

Supports flexible and ultra low crossing profile with thin distal shaft

* CrossFlex²™ technology available on 3.5 mm - 5.0 mm balloon sizes.

** Slim Seal technology available on 2.0 mm - 3.75 mm balloon sizes.

† Instead of hypotube with skive.

†† Similar to TREK™ RX.

Information contained herein for DISTRIBUTION outside of the U.S. ONLY.

Check the regulatory status of the device in areas where CE marking is not the regulation in force.

©2019 Abbott. All rights reserved. AP2947192-OUS Rev. A



TREK™ ORDERING INFORMATION

Offers a broad size matrix for more treatment options

		Balloon Length (mm)						
Balloon Diameter (mm)		6	8	12	15	20	25	30
TREK™ RX	2.25	1012271-06	1012271-08	1012271-12	1012271-15	1012271-20	1012271-25	1012271-30
	2.50	1012272-06	1012272-08	1012272-12	1012272-15	1012272-20	1012272-25	1012272-30
	2.75	1012273-06	1012273-08	1012273-12	1012273-15	1012273-20	1012273-25	1012273-30
	3.00	1012274-06	1012274-08	1012274-12	1012274-15	1012274-20	1012274-25	1012274-30
	3.25	1012275-06	1012275-08	1012275-12	1012275-15	1012275-20	1012275-25	1012275-30
	3.50	1012276-06	1012276-08	1012276-12	1012276-15	1012276-20	1012276-25	1012276-30
	3.75	1012277-06	1012277-08	1012277-12	1012277-15	1012277-20	1012277-25	1012277-30
	4.00	1012278-06	1012278-08	1012278-12	1012278-15	1012278-20	1012278-25	1012278-30
	4.50	-	-	1012279-12	1012279-15	-	-	-
	5.00	-	-	1012280-12	1012280-15	-	-	-
TREK™ OTW	2.25	1012404-06	1012404-08	1012404-12	1012404-15	1012404-20	1012404-25	1012404-30
	2.50	1012405-06	1012405-08	1012405-12	1012405-15	1012405-20	1012405-25	1012405-30
	2.75	1012406-06	1012406-08	1012406-12	1012406-15	1012406-20	1012406-25	1012406-30
	3.00	1012407-06	1012407-08	1012407-12	1012407-15	1012407-20	1012407-25	1012407-30
	3.25	1012408-06	1012408-08	1012408-12	1012408-15	1012408-20	1012408-25	1012408-30
	3.50	1012409-06	1012409-08	1012409-12	1012409-15	1012409-20	1012409-25	1012409-30
	3.75	1012410-06	1012410-08	1012410-12	1012410-15	1012410-20	1012410-25	1012410-30
	4.00	1012411-06	1012411-08	1012411-12	1012411-15	1012411-20	1012411-25	1012411-30
	4.50	-	-	1012412-12	1012412-15	-	-	-
	5.00	-	-	1012413-12	1012413-15	-	-	-

COMPLIANCE CHART

		Balloon Diameter (mm)									
Inflation Pressure (atm)		2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.50	5.00
2		1.87	2.15	2.16	2.49	2.91	2.93	3.06	3.24	3.69	4.10
3		1.95	2.20	2.35	2.56	2.98	3.04	3.22	3.46	3.91	4.28
4		2.00	2.27	2.41	2.66	3.05	3.11	3.31	3.56	4.04	4.43
5		2.06	2.34	2.49	2.73	3.12	3.20	3.41	3.67	4.17	4.59
6		2.12	2.40	2.56	2.81	3.19	3.28	3.50	3.78	4.29	4.73
7		2.17	2.45	2.63	2.87	3.25	3.38	3.59	3.89	4.40	4.86
8 NP		2.21	2.51	2.70	2.93	3.31	3.45	3.67	3.99	4.49	4.98
9		2.25	2.55	2.75	2.97	3.36	3.51	3.74	4.06	4.56	5.05
10		2.28	2.58	2.80	3.01	3.39	3.56	3.79	4.11	4.62	5.13
11		2.31	2.60	2.83	3.04	3.42	3.60	3.84	4.18	4.68	5.20
12		2.33	2.62	2.86	3.07	3.45	3.64	3.87	4.23	4.73	5.25
13		2.35	2.64	2.88	3.09	3.48	3.67	3.91	4.28	4.80	5.31
14 RBP		2.37	2.67	2.91	3.11	3.50	3.70	3.94	4.32	4.86	5.36
15		2.39	2.69	2.92	3.14	3.53	3.73	3.99	4.38	4.92	5.41
16		2.41	2.71	2.96	3.16	3.56	3.75	4.02	4.44	4.98	5.50
17		2.43	2.73	2.99	3.18	3.58	3.78	4.05	4.49	5.04	5.56

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for DISTRIBUTION outside of the U.S. only.

Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

www.CardioVascular.Abbott

©2019 Abbott. All rights reserved. AP2947192-OUS Rev. A

