

XIENCE Skypoint™

Everolimus Eluting Coronary Stent System

EXPANDED RANGE OF OTHER EXISTING DES



BETTER EXPANSION



EXCELLENT DELIVERABILITY




NEW 48 MM DESIGN



TRUSTED PATIENT OUTCOMES


SMALL STENT EXPANDS TO 3.75 MM



Test(s) performed and data on file at Abbott.

Small Stent Design

MEDIUM STENT EXPANDS TO 5.75 MM



Test(s) performed and data on file at Abbott.

Medium Stent Design

MULTIPLE PLATFORMS WITH POST EXPANSION RANGE

DES	SIZES	MAXIMAL EXPANSION
ONYX ^{(1) ‡}	2.0 - 2.5 mm	3.5 mm
	2.75 - 3.0 mm	4.0 mm
	3.5 - 4.0 mm	5.0 mm
	4.5 - 5.0 mm	6.0 mm
SYNERGY ^{(2) ‡}	2.25 - 2.75 mm	3.5 mm
	3.0 - 3.5 mm	4.25 mm
	4.0 - 5.0 mm	5.75 mm
ULTIMASTER TANSEI ^{(3) ‡}	2.25 - 3.0 mm	4.5 mm
	3.5 - 4.0 mm	5.5 mm
ORSIRO ^{(4) ‡}	2.25 - 3.0 mm	3.5 mm
	3.5 - 4.0 mm	4.5 mm

(1) RESOLUTE ONYX[®] BROCHURE - 02-2021 (2) SYNERGY[®] IFU 2019-06 (3) ULTIMASTER TANSEI[®] IFU 2020-03 (4) ORSIRO[®] IFU 2021-02

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at www.vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

Information contained herein for DISTRIBUTION in Europe, Middle East & Africa ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

Abbott Vascular International BVBA
Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.
‡ Indicates a third-party trademark, which is property of its respective owner.

www.cardiovascular.abbott

©2022 Abbott. All rights reserved. MAT-2202432 v1.0

